

# Ray·Ran

Part of the



# PRODUCT CATALOGUE

ISSUE 3

## CONTENTS

### RAW MATERIAL EVALUATION

#### MELT FLOW MEASUREMENT

|                           |   |
|---------------------------|---|
| RR/MFR100.....            | 3 |
| RR/MFR200.....            | 4 |
| RR/6MPCA.....             | 6 |
| TECHNI-TEST SOFTWARE..... | 7 |

#### MOISTURE ANALYSERS

|                 |    |
|-----------------|----|
| RR/MB23/25..... | 9  |
| RR/MB90.....    | 10 |
| RR/MB120.....   | 12 |

#### ANALYTICAL & PRECISION WEIGHING BALANCES

|                 |    |
|-----------------|----|
| PIONEER.....    | 13 |
| ADVENTURER..... | 14 |
| EXPLORER.....   | 15 |
| DISCOVERY.....  | 16 |

#### AUTO DENSITY GRADIENT APPARATUS

|                                   |    |
|-----------------------------------|----|
| RR/DGA (3 COLUMN & 6 COLUMN)..... | 18 |
| RR/DGA/CFS.....                   | 20 |
| TECHNI-TEST SOFTWARE.....         | 21 |

#### APPARENT BULK DENSITY

|  |    |
|--|----|
| ASTM D1895 METHOD A<br>& METHOD B.....       | 23 |
| ASTM D1895 METHOD C<br>& ISO METHOD R60..... | 24 |

#### DENSITY BALANCES

|                     |    |
|---------------------|----|
| ADVENTURER PRO..... | 27 |
| EXPLORER PRO.....   | 28 |
| DISCOVERY PRO.....  | 29 |

### TEST SAMPLE PREPARATION

#### CNC SAMPLE PROFILER

|                            |    |
|----------------------------|----|
| RR/CNC1.....               | 31 |
| RR/CNC2.....               | 32 |
| FLASHCUT CNC SOFTWARE..... | 33 |

#### TEST SAMPLE MOULDING APPARATUS

|              |    |
|--------------|----|
| RR/TSMA..... | 34 |
|--------------|----|

#### TEST SAMPLE CUTTING PRESSES

|             |    |
|-------------|----|
| RR/HCP..... | 37 |
| RR/PCP..... | 38 |

#### TEST SAMPLE CUTTERS

|             |    |
|-------------|----|
| RR/TSC..... | 40 |
|-------------|----|

#### AUTOCYCLE NOTCHING CUTTER

|            |    |
|------------|----|
| RR/NC..... | 42 |
|------------|----|

#### MICRO CUTTING MACHINE

|             |    |
|-------------|----|
| RR/MCM..... | 44 |
|-------------|----|

### SAMPLE AND COMPONENT TESTING

#### ADVANCED STATIC & DYNAMIC FRICTION TESTER

|                           |    |
|---------------------------|----|
| RR/FT.....                | 46 |
| TECHNI-TEST SOFTWARE..... | 48 |

#### FALLING DART IMPACT TESTER

|                           |    |
|---------------------------|----|
| RR/FDT-MB.....            | 50 |
| RR/FDT-MA.....            | 52 |
| RR/FDT-A1.....            | 54 |
| RR/FDT-A2.....            | 56 |
| TECHNI-TEST SOFTWARE..... | 57 |

#### PENDULUM IMPACT TESTER

|                           |    |
|---------------------------|----|
| RR/IMT.....               | 59 |
| TESTING METHODS.....      | 60 |
| TECHNI-TEST SOFTWARE..... | 61 |
| WITH PUNCTURE IMPACT..... | 62 |

#### FALLING WEIGHT IMPACT TESTER

|             |    |
|-------------|----|
| RR/FWT..... | 64 |
|-------------|----|

#### ADVANCED HDT/VICAT SOFTENING POINT

|                     |    |
|---------------------|----|
| RR/HDV2 & HDV4..... | 69 |
| RR/HDV6.....        | 72 |

#### COLD FLEX TESTER

|             |    |
|-------------|----|
| RR/CFT..... | 75 |
|-------------|----|

#### LOW TEMPERATURE BRITTLINESS TESTER

|             |    |
|-------------|----|
| RR/LTB..... | 77 |
|-------------|----|

#### ENVIRONMENTAL STRESS CRACKING APPARATUS

|             |    |
|-------------|----|
| RR/ESC..... | 79 |
|-------------|----|

#### UNRESTRAINED LINEAR THERMAL FILM SHRINKAGE LIQUID IMMERSION METHOD

|             |    |
|-------------|----|
| RR/FSL..... | 81 |
|-------------|----|

#### UNRESTRAINED LINEAR THERMAL FILM SHRINKAGE HOT PLATE METHOD

|              |    |
|--------------|----|
| RR/FSHP..... | 83 |
|--------------|----|

#### UNIVERSAL TESTING

|                        |    |
|------------------------|----|
| RR/UTM X250.....       | 84 |
| RR/UTM X350.....       | 86 |
| RR/UTM X500.....       | 88 |
| WINTEST™ ANALYSIS..... | 90 |

#### DIGITAL MICROMETERS

|              |    |
|--------------|----|
| RR/DM-T..... | 92 |
| RR/DM-B..... | 93 |

#### INDENTATION HARDNESS

|                         |    |
|-------------------------|----|
| RR/STD226 & STD227..... | 94 |
|-------------------------|----|


#### DUROMETER BENCH STAND

|             |    |
|-------------|----|
| RR/DBS..... | 95 |
|-------------|----|

### DEMONSTRATION & CALIBRATION

|                             |    |
|-----------------------------|----|
| MACHINE DEMONSTRATIONS..... | 96 |
|-----------------------------|----|

|                                      |    |
|--------------------------------------|----|
| SERVICING, CALIBRATION & REPAIR..... | 97 |
|--------------------------------------|----|



## RAW MATERIAL EVALUATION

### Basic Melt Flow System

# MODEL MFR100

The **MFR100** is the basic model offered within the Ray-Ran Range of Melt Flow Indexers and conforms to "method A" testing of International Test Standards. The manually operated test procedure is very simple to undertake. Molten Polymer is extruded through a closely controlled orifice (die) from the apparatus using pre set conditions of temperature and pressure produced by a dead weight system. The extruded polymer is cut off manually and is then weighed. Using the time interval to extrude the polymer its flow rate over 10 minutes can easily be determined, thus giving the **Melt Flow Index (MFI)** or **Melt Mass Flow Rate (MFR)** in **g/10 min**.

The MFR series of Melt Flow Indexers are supplied as standard with a replaceable hardened steel cylinder liner which can be simply swapped if a replacement is needed, or for Hastelloy (optional) if testing corrosive materials is a requirement. A Tungsten Carbide test die and hardened steel piston along with 2.16 kg test load and tooling ancillaries compliment the apparatus to get you testing straight out the box.

A 16th DIN temperature controller displaying set value and actual value, and a PT100 Platinum Resistance Thermocouple accurately control the barrel temperature to 0.1 deg C with a resolution of 0.1 deg C. Dual zone heating gives accurate temperature deviation along the length of the barrel to international test standards. An integrated timer is also fitted to accurately time the extrusion of the material in seconds to enable you to make cuts more accurately so the MFR can be calculated. All documentation is supplied including a product user manual and a fully traceable calibration certificate.



Optional weights up to and including 21.6kg can be supplied to cover all testing parameters to International Test standards and for the heavier weights the optional pneumatic weight loader can fitted to the apparatus for ease in loading the piston. The weight loader can also be used to hold a cleaning tool so cleaning can be done more simply. The improved die release feature makes removing the test die a simple task.

**Please be aware that this machine has no computer connectivity or file capture software functionality.**

Although the MFR value is not a fundamental property of the polymer, it does however, give an indication of the flow characteristics of the polymer and it has become one of the most widely used references for the quality control of polymers.

The machine is available in either 220-240v 50hz or 110v 60hz

## MFR100 BASIC MELT FLOW SYSTEM (RR/MFR100)

### TECHNICAL SPECIFICATION

- 16<sup>th</sup> DIN Digital temperature control
- Large Twin Colour Display
- Easy Set Point Operation
- Temperature Accurate to +/- 0.1°C
- Temperature Range 0 to 500°C
- Temperature Resolution +/- 0.1°C
- Digital Timer - Large easy to read display
- Count Range - 1 second to 99999 hours 59 seconds (8 character display)
- Tungsten carbide Test die, Hardened Steel Piston @ 0.325kg
- 2,16kg test load supplied as standard
- Simple replaceable cylinder liner design.
- Filling and cleaning tools included as standard
- Manual cutting
- Product user manual
- Traceable Calibration Certificate
- CE declaration certificate
- 1-year return to base warranty
- Available in 220-240V 50hz and 110-120V 60 hz, 10amp
- Conforms to Method A of International Test Standards ASTM D1238 & ISO1133 amongst others

### OPTIONAL ANCILLARIES:

- Pneumatic Weight Loader 150psi (11 bar) max
- Full Range of Test Weights Available from 1Kg to 21.6Kg
- Hastalloy cylinder, liner and die for corrosive materials

| WEIGHTS & DIMENSIONS: MODEL MFR100 |    |   |  |
|------------------------------------|----|---|--|
| Net Weight (kg)                    | 35 | Please note that shipping weights will increase if heavy test weights are included. |  |
| Width (cm)                         | 57 |   |  |
| Depth (cm)                         | 58 | For a complete set of test weights add approximately 30 kg to the weight.           |  |
| Height (cm)                        | 70 |   |  |

MFR200 BASIC MELT FLOW SYSTEM (RR/MFR200)

TECHNICAL SPECIFICATION

- Dedicated Microprocessor Control
- Dual Zone Heating
- Temperature Accuracy +/- 0.1°C
- Temperature Range ambient to 500°C
- Temperature Resolution +/- 0.01°C
- ECO mode
- Auto cut off device fitted with replaceable blades
- Intuitive 4" colour touch screen panel
- Tungsten Carbide Test Die, Hardened Steel Piston @ 0.325kg
- 2.16kg Test Weight supplied as standard
- Simple replaceable cylinder liner design
- Filling and cleaning tools included as standard
- Product user manual
- Traceable Calibration Certificate
- CE declaration certificate
- 1-year return to base warranty
- Available in 220-240V 50hz and 110-120V 60 hz, 10amp
- Conforms to Method A of International Test Standard ASTM D1238 & ISO1133 amongst others

TESTING FEATURES MFI:

- Temperature SV input
- Temperature AV display
- Number of cuts
- Time between cuts
- Mean or individual weight entry
- Auto MFI calculation
- Density @ melt temperature test and auto calculation

OPTIONAL ANCILLARIES

- Pneumatic Weight Loader 150psi (11 bar) max
- Full range of test weights available from 1kg to 21.6kg
- Hastelloy cylinder, liner and die for corrosive materials

| WEIGHTS & DIMENSIONS: MODEL MFR200 |    |   |
|------------------------------------|----|---|
| Net Weight (kg)                    | 35 | Please note that shipping weights will increase if heavy test weights are included. |
| Width (cm)                         | 57 |   |
| Depth (cm)                         | 58 | For a complete set of test weights add approximately 30 kg to the weight.           |
| Height (cm)                        | 70 |   |



RAW MATERIAL EVALUATION

Basic Melt Flow System

MODEL MFR200

The MFR200 melt flow indexer offers a more automated procedure to the MFR100 vastly improving results accuracy for “method A” testing requirements.

The test procedure is very simple to conduct. Molten Polymer is extruded through a closely controlled orifice (die) from the apparatus using pre set conditions of temperature and pressure produced by a dead weight system. The extruded polymer is cut off automatically and is then weighed. Using the time interval to extrude the polymer its flow rate over 10 minutes can easily be determined. Test results are given as the **Melt Flow Index (MFI)**. A **density at melt temperature test** can also be conducted.

The MFR series of Melt Flow Indexers are supplied as standard with a replaceable hardened steel cylinder liner which can be simply swapped if a replacement is needed or for Hastelloy (optional) if testing corrosive materials is a requirement. A Tungsten Carbide test die and hardened steel piston along with 2.16 kg test load and tooling ancillaries compliment the apparatus to get you testing straight out the box.

The apparatus features a rotary “Auto Cut Device” feature for consistent cutting accuracy of the extrudate eliminating the need for manual cutting.

For simplicity of use the apparatus is fitted with an intuitive 4” colour touch screen panel and microprocessor interface so test parameters can easily be set. Features include, temperature set point value (SV), actual temperature value (AV), number



of cuts, cut time, auto MFI calculation and “ECO” mode. Simply weigh the cut samples and enter the weights as a mean or individual cuts and let the microprocessor calculate the MFI value.

The integrated temperature controller and PT100 Platinum Resistance Thermocouple accurately control the barrel temperature to 0.1 deg C with a resolution of 0.01 deg C. Dual zone heating gives accurate temperature deviation along the length of the barrel to international test standards.

RAW MATERIAL EVALUATION

All documentation is supplied including a product user manual and a fully traceable calibration certificate.

Optional weights up to and including 21.6kg can be supplied to cover all testing parameters to International Test standards and for the heavier weights the optional pneumatic weight loader can fitted to the apparatus for ease in loading the piston. The weight loader can also be used to hold a cleaning tool so cleaning can be done more simply. The improved die release feature makes removing the test die a simple task.

**Please be aware that this machine has no computer connectivity or file capture software functionality.**

Although the MFR value is not a fundamental property of the polymer, it does however, give an indication of the flow characteristics of the polymer and it has become one of the most widely used references for the quality control of polymers.

The machine is available in either 220-240v 50hz or 110v 60hz



6MPCA ADVANCED MELT FLOW SYSTEM (RR/6MPCA)

TECHNICAL SPECIFICATION

- Microprocessor Temp Control
- Dual Zone Heating
- Temperature Accurate to +/- 0.1°C
- Temperature Range 0 to 400°C
- Temperature Resolution +/- 0.1°C
- Digital Encoder Accurate to +/- 0.02mm
- Multi slicing feature for accurate flow curve analysis
- Intuitive Menu Prompts
- Test Die, Piston & 2.16kg Test Weight supplied as standard
- Filling and cleaning tools included as standard
- Conforms to ASTM D1238, ASTM D3364, ISO 1133, DIN 53735 and others
- Electrical characteristics: 110v@60hz and 220v@50hz – fuse rating: 10amp
- **Techni-Test** software included as standard

TESTING FEATURES MFI:

- Operator list
- Material reference list
- Batch ref with data input
- Variable test temp input
- Variable test weight input
- Variable pre-heat input
- Material density data input
- Variable test distance
- Multi slicing feature
- High and low limit parameter setting

DENSITY:

- Operator list
- Material reference list
- Batch ref with data input
- Variable test temp input
- Variable test weight input
- Variable pre-heat input
- Variable test distance
- Numeric input of material weight in grams
- Automatic calculation of Density at test temperature

OPTIONAL ANCILLARIES:

- Weight Loader Option
- Full Range of Test Weights Available from 1Kg to 21.6Kg

WEIGHTS & DIMENSIONS:: MODEL 6MPCA

|                 |    |   |
|-----------------|----|---|
| Net Weight (kg) | 42 | Please note that shipping weights will increase if heavy test weights are included. |
| Width (cm)      | 57 |   |
| Depth (cm)      | 58 | For a complete set of test weights add approximately 30 kg to the weight.           |
| Height (cm)     | 70 |   |



Ray-Ran

RAW MATERIAL EVALUATION

Model 6 Advanced Microprocessor Controlled Melt Flow System with **Techni-Test** software

MODEL 6MPCA

The **6MPCA** is the most advanced model offered within the Ray-Ran Range of Melt Flow Indexers. The operating procedure is very simple to undertake using its on board advanced microprocessor technology. The large liquid crystal display (LCD) provides simple on-screen instructions reducing user error and the test parameters are easily entered via the alpha numeric membrane keypad. The apparatus can accurately determine results for **MFR**, **MVR** and **Density** at test temperature.

To conduct a test molten polymer is extruded through a closely controlled orifice (die) from the apparatus using user set test parameters of temperature and pressure produced by a dead weight system. Rotary encoder technology accurately determines piston displacement as the polymer is extruded and automatically starts the test when the piston is in the critical zone for testing.

Simple parameters are entered such as usernames, material reference numbers and batch numbers which are stored in lists for future recall and results presentation. Other parameters which are user defined are multi slicing and hi – lo limits. The multi slicing feature of the 6MPCA makes flow curve analysis instantly recognizable when the results are downloaded to the supplied **Techni-Test Software**. The operator inputs data for the number of slices to be made during the setup process



and the microprocessor accurately records the MFR result at each slice during the test. For quality control purposes the 6MPCA microprocessor can be set with High and Low limits which are clearly defined when the results are downloaded to a PC showing the user instantly if the material is a pass or fail.

If a materials density at test temperature is not known then a simple density test can be conducted by extruding molten

RAW MATERIAL EVALUATION

polymer over a selected test distance, weighing the cut off and inputting the weight into the microprocessor. The density result will be automatically calculated and stored for you to conduct your MFI test. After each test has been conducted the results are displayed on the LCD giving the MFI (g/10mins), Shear Stress (Pa), Shear Rate (1/sec), Viscosity (Pa/sec) and the Melt Volume Rate (cc/10min) which can then be downloaded via the on board mini USB connector or the Ethernet connection to Ray-Ran's dedicated **Techni-Test Software** where results are displayed in graphical and tabular form. CSV files of the results can also be saved which can be exported into other user programs for generating test reports.

The 6MPCA Melt Flow System is supplied as standard with a replaceable hardened steel cylinder liner, standard test die and piston along with 2.16 kg test load and cleaning ancillaries. The apparatus has dual Zone heating technology to keep the maximum allowable temperature variation along the length of the cylinder liner in accordance with the International test standard ISO1133

Optional weights can be supplied to cover all testing parameters to International Test standards and for the heavier weights the optional weight loader can fitted to the apparatus for ease in loading the piston. The weight loader can also be used as a hold back feature for materials which have a high flow rate. All documentation is supplied including **Techni-Test Software**, a product user manual and a fully traceable calibration certificate.

The machine is available in either 220-240v 50hz or 110v 60hz.



6MPCA ADVANCED MELT FLOW SYSTEM: TECHNI-TEST SOFTWARE

Ray-Ran

ENTER THE WORLD OF TECHNI-TEST

**Techni-Test** is an easy to use software package supplied with the **6MPCA Melt Flow system**, **Advanced Density Gradient Apparatus**, **Advanced Static & Dynamic Friction Tester**, **Advanced Falling Dart Impact Tester** and **Advanced Universal Pendulum Impact Tester**.

This allows the operator to monitor all aspects of the System during the test procedure.

Each tabular result displays the required information and results which are automatically calculated by the **Techni-Test Software** without the need for specialised equipment. Batch statistics such as Mean and Co-efficient of Variation (COV) are also displayed and are updated after each test result is downloaded. Test reports can be printed from the main screen when each test has been completed.

From the graph each multi-slice point is clearly identified giving accurate flow curve analysis of the material under test and by placing the cursor over each point the MFI value can be read on the screen. For Internal Quality Control Procedures at a glance, high and low limits are clearly displayed on the graph showing instantly if the material is a pass or fail. Multiple tests are clearly visible on the graph for results comparison within the batch and are highlighted in different colours for viewing.

In Test Results Viewer mode, users can upload saved results from previous tests for Material Comparison,



Data Manipulation or File Export. Abnormal results caused by air pockets within the molten material for example are clearly identified and can be removed from the test data bringing the batch statistics into a normal range ensuring that the test procedure does not have to be repeated saving time and material. Exporting the results file in Viewer Mode is simple. The Export file format is .CSV and can be opened with Microsoft Excel.





## RAW MATERIAL EVALUATION



# BASIC MOISTURE ANALYSERS MB23 AND MB25

Speed and Simplicity with Affordability

The **MB23** and **MB25 Basic Moisture Analysers** combine a high quality and durable ABS construction into a sleek, compact design and offer dependable, accurate results to quickly, effectively and affordably measure moisture content for a wide variety of applications.

Accurate results are provided on the custom-built backlit LCD screen at either 0.01g or 0.005g readability depending on the model and come with heating elements of either infrared in the MB23 or halogen in the MB25 for moisture determination between 50°C and 160°C. Each Moisture Analyser is ideal for textiles, wastewater, ceramics, food, and other applications.

Setup and operation are extremely easy. Simply press and hold the temperature or time buttons to set manual, automatic or timed duration drying parameters, add your sample and begin the test process. Results are displayed on the LCD at the end of the test as % moisture or % solids or weight (g), temperature and time and can be printed to an optional thermal printer via the onboard RS232 serial interface connection. This makes the MB23 and MB25 ideal for routine Quality Control tasks performed by operators of all skill levels.



Both the MB23 and MB25 Moisture Analysers have a compact footprint of only 17 x 13 x 18cm and have both been designed to take up less space and their easy to clean heating chamber is operator friendly for routine, inexpensive maintenance. Both the MB23 and MB25 come with a 90mm diameter sample pan and have a 110g weighing capacity.

## BASIC MOISTURE ANALYSER (RR/MB23/25)

### TECHNICAL SPECIFICATION

#### MB23

- MB23 Repeatability (Std Dev) (g) 0.3% for 3g sample/ 0.2% for 10g sample
- MB23 Readability 0.1% / 0.01g
- MB23 Infrared Heating

#### MB25

- MB25 Repeatability (Std Dev) (g) 0.2% for 3g sample/ 0.05% for 10g sample
- MB25 Readability 0.05% / 0.005g
- MB25 Halogen Heating

### STANDARD ON BOTH

- Weighing Capacity 110g
- Sample size 3g to 20g typical, 0.5g minimum
- RS232 Interface
- Timer 1 – 99 minutes, 30 second increments to 60 minutes
- Temperature range 50°C – 160°C in 5°C increments
- Electrical characteristics 100 – 240 VAC, 50/60 Hz
- LCD backlit display
- 50 extra sample pans included as standard
- CE, CSA/UL and FCC approved
- Sample Pan Handler (not supplied with MB23)

### OPTIONAL ANCILLARIES

- Security Locking Cable
- Security Lock
- Temperature Calibration Kit
- 50g M1 Calibration Weight
- Sample Pans (80/Box)
- Glass Fibre Pads (200/Box)
- Reusable Sample Pan (3/Box)
- Sample Pan Handler (MB23 only)
- Thermal Printer EU
- Thermal Printer UK
- Cable
- Software

### WEIGHTS & DIMENSIONS: MODEL MB25

|                 |     |
|-----------------|-----|
| Net Weight (kg) | 2.3 |
| Width (cm)      | 17  |
| Depth (cm)      | 28  |
| Height (cm)     | 13  |

TECHNICAL SPECIFICATION

MB90

- Capacity 90g
- Repeatability (Std Dev) (g) 0.18% for 3g sample/ 0.02% for 10g sample
- Readability 0.01% / 0.001g
- 40°C - 200°C
- Moisture range 0.01% - 100% (0.01% - 1000% for regain mode)
- Halogen heating
- Sample pan size 90mm
- RS232/USB Interface
- Operating Temperature range 10 - 40°C
- Electrical characteristics 100 – 240 VAC, 50/60 Hz
- QVGA, TFT Touch screen display
- Results in %moisture, %solids, %regain, time, temperature, weight, method name, drying curve and statistics
- CE, CSA/UL and FCC approved

OPTIONAL ANCILLARIES

- In use cover
- Scoop
- Temperature Calibration Kit
- Anti theft device
- Aluminium Sample Pans (80/Box)
- Filter Glass Fibre(200/Box)
- Reusable Sample Pan 7mm (3/Box)
- Reusable Sample Pan 14mm (3/Box)
- Impact Printer
- Thermal Printer UK

WEIGHTS & DIMENSIONS: MODEL MB90

|                 |     |
|-----------------|-----|
| Net Weight (kg) | 5.3 |
| Width (cm)      | 21  |
| Depth (cm)      | 36  |
| Height (cm)     | 18  |



Ray·Ran

RAW MATERIAL EVALUATION

# ADVANCED MOISTURE ANALYSERS MB90 AND MB120

Simple, Powerful, Affordable

The MB90 and MB120 Professional Moisture Analysers combine a high quality, sleek, compact design that effectively and affordably measure moisture content for a wide variety of applications.

Geared to increase precision and efficiency, Moisture Analysers feature a precisely controlled halogen heating system which helps to yield fast and repeatable results, with readability up to 0.01%. The robust housing and die-cast structure is designed for easy, tool-free cleaning and durability. You can access all of the MB120 and MB90's features and parameters on the easy-to-read colour touch screen display. The intuitive icon-based navigation guides users through every step of the operation process

The Temperature Guide helps to analyse a sample and determine the optimal drying temperature. A precisely controlled halogen heating system distributes heat uniformly over the sample, ensuring fast and repeatable results. Setup and operation are extremely easy. User selectable drying profiles for easy one-step drying and shut-off criteria to automatically end the test once moisture loss is no longer detected can be programmed. Test results



are displayed in %moisture, %solids, %regain, time, temperature, weight, method name, drying curve and statistics

Designed with convenient maintenance in mind, no tools are required for cleaning these Moisture Analysers. Simply remove the glass and tray to clean the sample chamber. The robust housing and die-cast base guarantee durability

TESTING MADE SIMPLE

The MB Series features a wide, bright touch screen display and intuitive menu navigation for fast, accurate moisture analysis.

- **Touch screen Display:** Easy to read and understand
- **Icon-Driven Menu Navigation:** Easy to follow, with no training required
- **4 Mechanical Keys for Fast One-Step Operation:** Tare, Start/Stop, Print, On/Off
- **USB and RS232 Ports:** Quick connections to a PC or flash drive



Ray·Ran

RAW MATERIAL EVALUATION

MOISTURE ANALYSIS  
MADE SIMPLE: JUST 3  
EASY STEPS FOR ROUTINE  
SAMPLE TESTING!

1. **Tare empty sample pan**
2. **Add sample according to best practices**
3. **Press start and see your results!**

EXCEPTIONAL PRECISION

The MB120 and MB90 provide precise results for more confident testing.

- **Fast Infrared Halogen Heating Technology:** Halogen heating dries samples in seconds and performs up to 40% faster than traditional Infrared methods.
- **Uniform Drying for Precise, Repeatable Results:** The reflective sample chamber ensures uniform heating of each sample to yield reliable results of 0.01% readability.
- **Flexible and Powerful:** Features a full operating temperature range of 40–200 °C (MB90) or 40–230 °C (MB120) in 1-degree increments.
- **Four Drying Profiles to Ensure Accuracy (MB120)** Customize heating of specific samples with temperature control for one-step drying—set temperature and time and begin!  
**FAST:** Quickest possible time to temperature  
**STANDARD:** Minimal overshoot while achieving temperatures quickly  
**RAMP:** Controlled ramping to temperatures



**STEP:** Set three temperatures at three time intervals

- **Shut-Off Criteria for Automatic Test Completion (MB120):** Customize your own test endpoint criteria to automatically end the test once moisture loss is no longer detected. Choose from up to 7 criteria—Manual, Timed, Auto Free percent weight loss per second, Auto Free mg loss per

ADDITIONAL FEATURES—MB90

ABS top housing, stainless steel pan support, stainless steel pan handler, in-use cover, 3 shut-off criteria (manual, timed or auto), 2 heating profiles (standard and fast), storage up to 2 methods and user guide, 40–200 °C heating temperature range (1 °C increments), 13 operational languages (display only, no key in)

ADDITIONAL  
FEATURES—MB120

ABS top housing, stainless steel pan support, stainless steel metal pan handler, in-use cover, 7 shut-off criteria (manual, timed, auto free %/s, auto free mg/s, A30, A60 and A90), 4 heating profiles (standard, fast, ramp and step), user guide, 40–230 °C heating temperature range (1 °C increments), 13 operational languages (display only, no key in)

TECHNICAL SPECIFICATION

MB120

- Capacity 120g
- Repeatability (Std Dev) (g) 0.05% for 3g sample/ 0.015% for 10g sample
- Readability 0.01% / 0.001g
- Temperature range 40°C - 230°C
- Moisture range 0.01% - 100% (0.01% - 1000% for regain mode)
- Halogen heating
- Sample pan size 90mm
- RS232/USB Interface
- Operating Temperature range 10 - 40°C
- Electrical characteristics 100 – 240 VAC, 50/60 Hz
- QVGA, TFT Touch screen display
- Results in %moisture, %solids, %regain, time, temperature, weight, method name, drying curve and statistics
- CE, CSA/UL and FCC approved

OPTIONAL ANCILLARIES

- In use cover
- Scoop
- Temperature Calibration Kit
- Anti theft device
- Aluminium Sample Pans (80/Box)
- Filter Glass Fibre(200/Box)
- Reusable Sample Pan 7mm (3/Box)
- Reusable Sample Pan 14mm (3/Box)
- Impact Printer
- Thermal Printer UK

WEIGHTS & DIMENSIONS: MODEL MB120

|                 |     |
|-----------------|-----|
| Net Weight (kg) | 5.3 |
| Width (cm)      | 21  |
| Depth (cm)      | 36  |
| Height (cm)     | 18  |



RAW MATERIAL EVALUATION

PIONEER  
WEIGHING BALANCE

The **Pioneer** series of analytical and precision balances are designed for basic routine weighing in a variety of laboratory, industrial and education applications.

With the right combination of performance and features the Pioneer offers uncomplicated performance for all your basic weighing needs and includes parts counting and percentage weighing.

The integral up-front level indicator quickly ensures the

balance is level prior to use and to help adapt the balance to real life working conditions the Pioneer is designed with selectable environmental settings. When working in harsh environments, the balance can be adjusted to compensate for vibrations and other disturbances. Conversely, the balance can be adjusted for working in slow filling applications where sensitivity is a must.

WEIGHING CAPACITIES

- Analytical** – 65g/0.1mg up to 210g/0.1mg available
- Precision** – 210g/1mg up to 4100g/0.1g available

OTHER STANDARD FEATURES  
AND EQUIPMENT

AC adapter, user selectable span calibration points, software reset menu, stability indicator, Auto Tare, user selectable communications settings, user selectable printing options and stainless-steel platform.

The **Pioneer** can be supplied with OIML approval.



PIONEER: WEIGHING BALANCE

TECHNICAL SPECIFICATION

- Weighing Capacity 65g up to 4100g
- Repeatability (g) 0.0001 up to 0.1
- Readability (g) 0.0001 up to 0.1
- Liquid Crystal Display (LCD )
- RS232 with GLP/GMP Data Output
- Rigid Metal/ABS construction
- Glass draft-shield with flip top door
- AC adapter included
- Multiple Application Modes
- Selectable Environmental Settings
- 18 selectable units of measurement (non OIML approved)

OPTIONAL ANCILLARIES

- Thermal Printer
- Impact Printer
- Printer Cable
- Computer software
- Display cover
- Internal calibration
- Traceable Calibration Certificate

| WEIGHTS & DIMENSIONS: PIONEER |     |
|-------------------------------|-----|
| Net Weight (kg)               | 4.5 |
| Width (cm)                    | 20  |
| Depth (cm)                    | 32  |
| Height (cm)                   | 9.2 |

TECHNICAL SPECIFICATION

- Weighing Capacity 110g up to 8100g
- Repeatability (g) 0.0001 up to 0.1
- Readability (g) 0.0001 up to 0.1
- 2-Line Liquid Crystal Display (LCD ) with Backlight
- RS232 with GLP/GMP Data Output
- Rigid Metal/ABS construction
- Glass draft-shield with flip top door
- AC adapter included
- Multiple Application Modes
- Selectable Environmental Settings
- 18 selectable units of measurement (non OIML approved)

OPTIONAL ANCILLARIES

- Thermal Printer
- Impact Printer
- Printer Cable
- Computer software
- Display cover
- Internal calibration
- Traceable Calibration Certificate

| WEIGHTS & DIMENSIONS: ADVENTURER |     |
|----------------------------------|-----|
| Net Weight (kg)                  | 4.4 |
| Width (cm)                       | 22  |
| Depth (cm)                       | 30  |
| Height (cm)                      | 7.2 |

ADVENTURER WEIGHING BALANCE

The **Adventurer** is one of the industry’s most versatile balances – it counts, it sums, it holds and it delivers traceable results.

With a full range of **Analytical** and **Precision models**, no other balance in its class offers all these features in a compact design. Internal calibration ensures weighing accuracy by allowing you to calibrate the balance prior to use without the need of external weights and the environmental settings offer three filter modes and adjustable zero tracking which are ideal for working in slow filling applications where sensitivity is a must. When working in harsh environments, the user



14

can adjust balance settings to compensate for vibrations and other disturbances. All Adventurer models come with a **built-in RS-232 serial interface**, and with **full GLP/GMP data output**, it is sure to meet all your traceability and compliance requirements.

OTHER FEATURES INCLUDE

- Parts Counting with Optimisation** – Allows the user to count based on a calculated piece weight. With optimisation, counting accuracy is improved through auto re-calculation of the piece weight
- Percentage Mode** – Allows the user to display the weight of a sample as a percentage of a reference weight
- Check weighing Mode** – Allows the user to check the weight of a sample against a pre-set target
- Totalization Mode** – Allows the user to sum a series of sample weights
- Display Hold Mode** – Allows the user to hold the highest stable weight on the display until cleared

WEIGHING CAPACITIES

- Analytical** – 110g/0.1mg up to 260g/0.1mg available
- Precision** – 210g/1mg up to 8100g/0.1g available

**The Adventurer can be supplied with OIML approval.**

EXPLORER WEIGHING BALANCE

The **Explorer** series of analytical and precision balances has established itself as market leader in its class.

A compact design, the Explorer delivers accurate results within seconds improving operator efficiency, productivity and throughput, with a stabilisation time up to **50% faster than basic models**. Its superior anti-vibration system makes the Explorer **ideal for unstable working conditions**. It is supplied with a frameless, flip top anti-static glass draft-shield which provides unobstructed access to the working chamber. Simple prompts using Icons on the **large colour VGA display** make navigating the menus easy and the resistive touch screen display quickly responds to the operators touch or stylus.



Also included on the Explorer is a four-point touch-less sensor for total hands-free operation and the internal calibration procedure ensures accuracy is maintained even if ambient room temperatures effect the weighing accuracy. Test results can be downloaded via the onboard RS232 serial interface connector direct to Excel or to an optional thermal printer.

OTHER FEATURES INCLUDE

Qwerty keyboard and numeric keypad, 14 built in applications, 11 operating languages, NTEP and OIML models available,

menu lock switch, security bracket, integral weigh below hook for below balance weighing applications, removable stainless steel weighing platform, stability indicator, overload and under load indicators, auto standby Weighing (17 units + custom units), Parts Counting, % weighing, Animal/Dynamic Weighing, Check weighing, Filling, Totalization, Formulation, Differential weighing, Density determination, Peak hold, Ingredient costing, Pipette adjustment, Gross/Net/Tare Weighing

**The Explorer can be supplied with OIML approval.**

TECHNICAL SPECIFICATION

- Weighing Capacity 120g up to 10200g
- Repeatability (g) 0.0001 up to 0.1
- Readability (g) 0.0001 up to 0.1
- Auto Calibration
- RS232 Interface & USB
- Display Size 5.7 in / 145 mm (diagonal)
- Display Full-colour VGA graphic display, 4-wire resistive touch screen
- Electrical characteristics AC Adapter Input: 100-240 VAC 0.6A 50-60 Hz
- CE, CSA/UL and FCC approved
- 18 selectable units of measurement (non OIML approved)

OPTIONAL ANCILLARIES

- Thermal Printer
- Impact Printer
- Printer Cable
- Ethernet Interface Connector

| WEIGHTS & DIMENSIONS: EXPLORER |     |
|--------------------------------|-----|
| Net Weight (kg)                | 6.9 |
| Width (cm)                     | 23  |
| Depth (cm)                     | 39  |
| Height (cm)                    | 36  |

15

TECHNICAL SPECIFICATION

- Analytical Weighing Capacity up to 110g up to 310g
- Analytical Repeatability (g) 0.0001
- Analytical Readability (g) 0.0001
- Semi Micro Weighing Capacity up to 81g up to 210g
- Semi Micro Repeatability (g) 0.00001 – 0.0001
- Semi Micro Readability (g) 0.00001 – 0.0001
- Auto Calibration
- RS232 Interface with full GLP/GMP Protocol
- Display Size 4 x 1 / 10 x 2.5
- 2-line Alphanumeric Backlit LCD Display
- External Adapter, 100-120VAC 150mA, 220-240VAC 100mA, 50/60Hz
- CE, CSA/UL and FCC approved
- 18 selectable units of measurement (non OIML approved)

OPTIONAL ANCILLARIES

- Thermal Printer
- Impact Printer
- Printer Cable
- Ethernet Interface Connector
- Computer Software
- Dust Cover
- Anti Theft Device

| WEIGHTS & DIMENSIONS: DISCOVERY |      |
|---------------------------------|------|
| Net Weight (kg)                 | 10.2 |
| Width (cm)                      | 20   |
| Depth (cm)                      | 46   |
| Height (cm)                     | 30   |

DISCOVERY  
WEIGHING BALANCE

The **Discovery** weighing Balance is the professional’s choice when it comes to weight measurement.

It features a rigid construction and an advanced internal calibration system offering the performance you would expect to find only in premium-priced balances. A compact design, the Discovery Pro delivers accurate results within seconds improving operator efficiency and productivity. The easy to use balance software utilizes text prompts to guide users through applications and balance setup. Delivered through a 2-line backlit LCD display and clearly marked buttons, the Discovery Pro tells you what is going on and what you must do next. It is supplied with a frameless, flip top anti-static glass draft-shield which provides unobstructed access to the working chamber. To minimize weighing errors and ensure accurate measurements, two internal weights are used to perform the internal linear calibration procedure which automatically calibrates the balance when it senses a temperature change significant enough to affect the weighing accuracy. The balance features a rear mount load cell for ultimate, accuracy and isolation.



The Discovery’s repeatability and linearity performance is 2x better than most semi-micro balances in its class. Made from a glass and steel construction it resists the effects of static and corrosion increasing balance stability and accuracy, as well as extending product life.

OTHER FEATURES INCLUDE

Selectable environmental filters, protective in-use cover, integral weigh below hook, easy to use keypad, up-front level indicator, stability indicator, mechanical and software overload/under load protection, AC adapter, user selectable span calibration points, auto tare, user selectable printing options, user selectable communications settings, user selectable data print options, user definable project and user ID’s, software reset menu, software lockout menu

WEIGHING FEATURES

Weighing, Parts Counting with Automatic Sample Recalibration, Percent Weighing, Check weighing, Animal/Dynamic Weighing, Gross/Net/Tare Weighing, Totalization, High Point, Density, Statistics and Pipette Calibration

**The Discovery can be supplied with OIML approval and is offered in analytical and semi-micro versions.**





### 3 COLUMN DGA & 6 COLUMN DGA (RR/DGA)

# AUTO DENSITY GRADIENT APPARATUS

The **Advanced Auto Density Gradient Apparatus** from Ray-Ran has become the world's benchmark for accurate density measurement of small solid specimens using the column method. Built with operator simplicity in mind, its ease of operation and high accuracy makes it ideal for product development and quality control within production, research and development labs and teaching institutions.

Offered as a **3 or 6 column version**, the built-in on-board Microprocessor System accurately calculates the specimen's density to a resolution of 0.0001g/ml more quickly and more accurately by measuring the samples position in the column relative to the calibrated marker floats. Once the correct position of the sample is recorded the density is displayed on the LCD screen. Column linearity can also be checked and shown in **Techni-Test**.

The on-board liquid crystal display (LCD) provides simple on-screen instruction and using the alpha/numeric keypad test parameters such as operator list, material reference numbers and batch numbers can be entered and stored for results presentation. Batch statistics of mean, standard deviation and co-efficient of variation are automatically updated after each test and results are shown in graphical and

## RAW MATERIAL EVALUATION

tabular format for analysis when downloaded to the supplied **Techni-Test Software**.

For results accuracy the Density Gradient Apparatus is supplied with built in digital temperature control to ensure the temperature of each column is at 23°C +/- 0.1°C. For ambient temperatures in excess of 23°C the integrated cooling coil can be used in conjunction with an optional water chiller unit so the correct temperature can be maintained. This can be even more accurately controlled using the optional microprocessor-controlled water flow valve kit.

Test results are displayed in g/cm<sup>3</sup> which can then be downloaded via the on-board USB or Ethernet connection to Ray-Ran's dedicated **Techni-Test Software** where results are displayed in graphical and tabular form. CSV files of the results can also be saved which can be exported into other user programs for generating test reports.

To build a column a peristaltic pumped filling system is supplied with the apparatus. Column filling speeds are variable from 0.5ltrs/hr to 1.5ltrs/hr. The resulting solution has a density which when introduced into the column increases uniformly from top to bottom. Calibrated glass marker floats with known densities are introduced into the column and sink to a point where

their density matches the solution. A linear encoder which is attached to a trammel guide is used to calibrate the gradient of the column by focusing the optical microscope onto the centre of the calibrated float and entering the floats density value into the microprocessor. Once test samples have been introduced into the column and have reached a point of static equilibrium they are accurately sighted using the microscope. The resultant density of the sample is obtained and displayed on the LCD without the need for complicated graphs. An automatic sweep mechanism is used to extract the floats and samples from each column.



### 3 COLUMN DGA & 6 COLUMN DGA (RR/DGA)

#### TECHNICAL SPECIFICATION

- Automatic density calculation
- Automatic calibration system
- LCD display
- Resolution 0.000001 g/ml
- Accuracy 0.0001 g/ml
- Backlight
- Variable speed pumped filling system
- Twin conical filling flasks
- Automatic magnetic stirrer
- Automatic sweep mechanism
- Stainless steel sweep baskets
- 7x optical microscope
- Digital temperature control to 0.1°C
- Mini USB port
- Cooling coil
- 110v 60hz and 240v 50Hz
- Product user manual
- Traceable calibration certificate
- CE declaration certificate
- 1 year return to base warranty
- Conforms to ISO1183 & ASTM D1505

#### OPTIONAL ANCILLARIES

- Thermal Printer
- Water chiller unit
- Glass marker floats 0.7g/ml up to 2.2g/ml
- Microprocessor controlled filling system

| WEIGHTS & DIMENSIONS:<br>AUTO DENSITY GRADIENT APPARATUS |              |              |
|--|--------------|--------------|
|  | 3 COLUMN DGA | 6 COLUMN DGA |
| Net Weight (kg)  | 70           | 80           |
| Width (cm)   | 85           | 120          |
| Depth (cm)   | 40           | 40           |
| Height (cm)  | 110          | 110          |

TECHNICAL SPECIFICATION

- Small footprint filler unit with basic operating panel
- LCD display
- Fully programmable column volume, and upper and lower densities
- More accurate and consistent than other filling methods
- Variable speed pumped filling system
- No need to premix liquids
- g/ml or kg/m3 units
- Variable top and bottom buffer volumes
- Previous stored filling routines
- Automatic density calculation
- Mains lead
- 5mtrs of 4mm tubing supplied
- Weighted end for the outlet tube
- Spare pump tubes & connectors
- 110v - 240v 50 – 60Hz
- Product user manual
- CE declaration certificate
- 1 year return to base warranty
- Product user manual

| WEIGHTS & DIMENSIONS: RR/DGA/CFS |     |
|----------------------------------|-----|
| Net Weight (kg)                  | 5.3 |
| Width (cm)                       | 32  |
| Depth (cm)                       | 25  |
| Height (cm)                      | 11  |

MICROPROCESSOR CONTROLLED DENSITY GRADIENT FILLING SYSTEM

The **Microprocessor Controlled Filling System** offered by Ray-Ran is by far the most accurate way of producing a density gradient within your column. Not only does it effortlessly integrate with the **Ray-Ran 3 and 6 column equipment** but can also be used with any other make of density gradient column. It saves considerable time in the preparation



and building of your density gradient and overcomes a lot of the trial and error that is associated with this procedure.

Buffer zones can be added to the top and bottom of the column to slow down the deterioration of the column and gradient settings can be stored for future use including controlled top and bottom densities, column volume, tube volume and buffer volumes as well as selectable density units, with either g/ml or kg/m<sup>3</sup> available.



The filling procedure is fully automated using indexing stepper motors and peristaltic pumps and any two miscible liquids can be used to give the heavy and light densities without the need to be pre-mixed prior to being pumped into the column.

The filling system is a faster, more accurate and more consistent method of building a density gradient column increasing time and productivity.

DENSITY MEASUREMENT MADE SIMPLE WITH TECHNI-TEST

**Techni-Test** is an easy to use software package supplied with the **Advanced Density Gradient Apparatus** which allows user defined test data to automatically download from the apparatus for report information and for the operator to analyse all aspects of the Density test.

From the graph each density reading is clearly identified giving accurate data point analysis of the materials under test. Placing the cursor over each plotted point the Density value can be read on the screen. Density results are also shown in tabular form. Test reports can be printed from the main screen when required by selecting the Print Report dialogue box to print the report. Batch statistics such as Mean, Standard Deviation and Co-efficient of Variation (COV)



are also displayed and are updated after each test result is downloaded. By Selecting the tabs you can select between the test results data and look at the linearity of the column that's been calibrated.

In Test Results Viewer mode users have the ability to upload saved results from previous test for Material Comparison, Data Manipulation or File Export. Abnormal results that could be caused by air or dust particles are clearly identified and can be removed from the test data bringing the batch statistics into a normal range ensuring that the test procedure does not have to be repeated saving time and material. Exporting the results file in Viewer Mode is simple. The Export file format is .CSV and can be opened with Microsoft Excel.



RAW MATERIAL EVALUATION

APPARENT BULK DENSITY

**Bulk density** is a property of powders, granules and other “divided” solids, especially used in reference to mineral components (soil, gravel), chemical substances, (pharmaceutical) ingredients, foodstuff or any other masses of corpuscular or particulate matter. It is defined as the mass of many particles of the material divided by the total volume they occupy. The total volume includes particle volume, inter-particle void volume and internal pore volume.

Bulk density is not an intrinsic property of a material; it can change depending on how the material is handled. For example, a powder poured into a cylinder will have a particular bulk density; if the cylinder is disturbed, the powder particles will move and usually settle closer together, resulting in a higher bulk density. For this reason, the bulk density of powders is usually reported both as “freely settled” (or “poured” density) and “tapped” density (where the tapped density refers to the bulk density of the powder after a specified compaction process, usually involving vibration of the container.)

ASTM D1895 METHOD A

Primarily used to measure the apparent density of fine granules that can be poured readily through a “V” shaped funnel, the material under test is allowed to flow into a cylindrical cup with a known volume of 100cm<sup>3</sup>.

TESTS THAT CAN BE CONDUCTED ARE:

- Apparent Density
- Bulk Factor
- Pourability



ASTM D1895 METHOD B

The larger of the ASTM family Method B is primarily used to measure the apparent density of larger coarse granular materials, dice or pellets that cannot be poured readily through the method A funnel. The material under test is allowed to flow into a cylindrical cup with a known volume of 400cm<sup>3</sup>.

TESTS THAT CAN BE CONDUCTED ARE:

- Apparent Density
- Bulk Factor
- Pourability



ASTM D1895 METHOD A & METHOD B

ASTM D1895 METHOD A TECHNICAL SPECIFICATION

- Conforms to ASTM D1895 method A
- V shaped funnel
- Measuring cup 100cm<sup>3</sup>
- Stand with funnel shut off

WEIGHTS & DIMENSIONS: ASTM D1895 A

|                 |    |
|-----------------|----|
| Net Weight (kg) | 4  |
| Width (cm)      | 18 |
| Depth (cm)      | 23 |
| Height (cm)     | 25 |

ASTM D1895 METHOD B TECHNICAL SPECIFICATION

- Conforms to ASTM D1895 method B
- V shaped funnel
- Measuring cup 400cm<sup>3</sup>
- Stand with funnel shut off

WEIGHTS & DIMENSIONS: ASTM D1895 A

|                 |    |
|-----------------|----|
| Net Weight (kg) | 9  |
| Width (cm)      | 24 |
| Depth (cm)      | 33 |
| Height (cm)     | 52 |

ASTM D1895 METHOD C & ISO METHOD R60

ASTM D1895 METHOD C  
TECHNICAL SPECIFICATION

- Conforms to ASTM D1895 method C
- Measuring funnel
- Weight plunger
- Lead shot

| WEIGHTS & DIMENSIONS: ASTM D1895 C |          |
|------------------------------------|----------|
| Net Weight (kg)                    | 4        |
| Width (cm)                         | 10       |
| Depth (cm)                         | 10       |
| Height (cm)                        | 30 (max) |

ISO METHOD R60  
TECHNICAL SPECIFICATION

- Conforms to ISO R60 specification
- Measuring cup 100cm<sup>3</sup>
- Stand with funnel shut off

| WEIGHTS & DIMENSIONS: ISO METHOD R60 |    |
|--------------------------------------|----|
| Net Weight (kg)                      | 4  |
| Width (cm)                           | 18 |
| Depth (cm)                           | 23 |
| Height (cm)                          | 25 |

ASTM D1895 METHOD C

For course flakes, strands, chips and cut fibres that cannot be poured using test methods A and B Ray-Ran offer a measuring cylinder and plunger to method C of the ASTM test standard. Since these types of materials to be tested are very bulky when loosely poured and are usually compressed to lessen the bulk, a measure of their density under a small compression load is very useful. For this test a measuring cylinder of 1000 cm<sup>3</sup> is supplied along with a scaled weight plunger with 1mm graduations on the outside. Lead shot is used to increase the plunger weight to 2300g to compress the material under test.

TESTS THAT CAN BE CONDUCTED ARE

- Apparent Density



ISO METHOD R60

Primarily used to measure the apparent density of moulding material that can be poured readily through a specified funnel, the material under test is allowed to flow into a cylindrical cup with a known volume of 100cm<sup>3</sup>.

TESTS THAT CAN BE CONDUCTED ARE

- Apparent Density
- Bulk Factor
- Pourability



OTHER  
MODELS ARE  
AVAILABLE ON  
REQUEST

PLEASE CONTACT  
RAY-RAN FOR MORE  
INFORMATION  
REGARDING YOUR  
SPECIFIC REQUIREMENT



RAW MATERIAL EVALUATION

# ADVENTURER PRO DENSITY DETERMINATION BALANCE

The **Adventurer Pro Analytical Series Density Determination Balance** has established itself as one of the industry’s most versatile balances. A compact design, it is supplied with a metal base and ABS shell, ABS top housing, glass draft-shield with flip top door and replaceable in-user cover. It determines the density of fluids and the specific gravity of solid materials by using the **Archimedeian Buoyancy Method**. Simple menu prompts with large backlit 2-Line liquid crystal display makes the equipment user friendly with a display resolution at 0.1mg. The internal calibration system ensures weighing accuracy allowing you to calibrate the balance prior to use without the need of external weights.

To conduct a test the sample is first weighed in air and again when submerged in a liquid. By determining the loss in weight which is

directly proportional to the mass of the fluid it has displaced, the specific gravity and density is readily calculated. Test results can be downloaded via the onboard RS232 serial interface connector to optional dedicated PC software or thermal printer.

OTHER FEATURES INCLUDE

User selectable environmental settings, stability indicator, mechanical and software overload/underload protection, AC adapter, user selectable span calibration points, auto tare, user selectable printing options, user selectable communications settings, user selectable data print options, user definable project and user ID’s, software reset menu, software lockout menu and integral weigh-below hook.

The Adventurer Pro is an essential, cost effective tool for Laboratory and Quality Control environments.

ADVENTURER PRO DENSITY BALANCE

TECHNICAL SPECIFICATION

- Weighing Capacity up to 260g
- Repeatability (g) 0.0001
- Readability (g) 0.0001
- 2-Line Liquid Crystal Display (LCD ) with Backlight
- RS232 with GLP/GMP Data Ouput
- Rigid Metal/ABS construction
- Glass draft-shield with flip top door
- AC adapter included
- Multiple Application Modes
- Selectable Environmental Settings

OPTIONAL ANCILLARIES

- Thermal Printer
- Impact Printer
- Printer Cable
- Computer software
- Internal Self Calibration

WEIGHTS & DIMENSIONS: ADVENTURER PRO

|                 |      |
|-----------------|------|
| Net Weight (kg) | 4.4  |
| Width (cm)      | 22   |
| Depth (cm)      | 30   |
| Height (cm)     | 30.5 |

TECHNICAL SPECIFICATION

- Weighing Capacity up to 320g
- Repeatability (g) 0.0001
- Readability (g) 0.0001
- Auto Calibration
- RS232 Interface & USB
- Display Size 5.7 in / 145 mm (diagonal)
- Display Full-colour VGA graphic display, 4-wire resistive touch screen
- Electrical characteristics AC Adapter  
Input: 100-240 VAC 0.6A 50-60 Hz
- CE, CSA/UL and FCC approved

OPTIONAL ANCILLARIES

- Thermal Printer
- Impact Printer
- Printer Cable
- Ethernet Interface Connector

| WEIGHTS & DIMENSIONS: EXPLORER PRO |     |
|------------------------------------|-----|
| Net Weight (kg)                    | 6.9 |
| Width (cm)                         | 23  |
| Depth (cm)                         | 39  |
| Height (cm)                        | 36  |

EXPLORER PRO DENSITY DETERMINATION BALANCE

The **Explorer Pro Analytical Series Density Determination Balance** has established itself a market leader in its class. A compact design, the Explorer Pro delivers accurate results within seconds improving operator efficiency, productivity and throughput, with a stabilisation time up to 50% faster than basic models.

Its superior anti-vibration system makes the Explorer Pro ideal for unstable working conditions. It is supplied with a frameless, flip top anti-static glass draft-shield which provides unobstructed access to the working chamber. Simple prompts using Icons on the large colour VGA display makes navigating the menus easy and the resistive touch screen display quickly responds to the operators touch or stylus. Also included on the Explorer Pro is a four-point touchless sensor for total hands-free operation and the internal calibration procedure ensure accuracy to 0.1mg is maintained even if ambient room temperatures affect the weighing accuracy.



The Explorer Pro accurately determines the density of fluids and the specific gravity of solid materials by using the Archimedean Buoyancy Method. To conduct a test the sample is first weighed in air and again when submerged in a liquid. By determining the loss in weight which is directly proportional to the mass of the fluid it has displaced, the specific gravity and density is readily calculated. Test results can be downloaded via the onboard RS232 serial interface connector direct to Excel or to an optional thermal printer.

OTHER FEATURES INCLUDE

Qwerty keyboard and numeric keypad, 14 built-in applications, fully-automatic internal calibration system, fast stabilization time, four touchless sensors, up to 11 operating languages, NTEP and OIML models available, menu lock switch, security bracket, integral weigh below hook for below balance weighing applications, removable stainless steel weighing platform, stability indicator, overload and underload indicators, auto standby.


The Explorer Pro is an essential, cost effective tool for Laboratory and Quality Control environments.

DISCOVERY PRO DENSITY DETERMINATION BALANCE

The **Discovery Pro Density Determination Balance**

is the professional's choice when it comes to measuring density. The Discovery Pro features a rigid construction and an advanced internal calibration system offering the performance you would expect to find only in premium-priced balances.

A compact design, the Discovery Pro delivers accurate results within seconds improving operator efficiency and productivity. The easy to use balance software utilizes text prompts to guide users through applications and balance setup. Delivered through a 2-line backlit LCD display and clearly marked buttons, the Discovery Pro tells you what is going on and what you must do next. It is supplied with a frameless, flip top anti-static glass draft-shield which provides unobstructed access to the working chamber. To minimize weighing errors and ensure accurate measurements, two internal weights are used to perform a linear calibration. The internal calibration system automatically calibrates the balance when it senses a temperature change significant enough to affect the weighing accuracy.



The Discovery Pro provides users with four methods of density determination to choose from depending on their need. They are:

1. Determine density for solids denser than water
2. Determine density for solids less dense than water
3. Determine density for liquid density (sinker needed – not included)
4. Determine density for Porous material (impregnated with oil).

After entering a couple of parameters into the balance, the built-in density table and mass measurement of the balance take over to provide rapid and accurate density results, with no manual calculations needed.

OTHER FEATURES INCLUDE

Selectable environmental filters, protective in-use cover, integral weigh below hook, easy to use keypad, up-front level indicator, stability indicator, mechanical and software overload/underload protection, AC adapter, user selectable span calibration points, auto tare, user selectable printing options, user selectable communications settings, user selectable data print options, user definable project and user ID's, software reset menu and software lockout menu.

TECHNICAL SPECIFICATION

- Weighing Capacity up to 310g
- Repeatability (g) 0.00001
- Readability (g) 0.00001
- Auto Calibration
- RS232 Interface with full GLP/GMP Protocol
- Display Size 4 x 1 / 10 x 2.5
- 2-line Alphanumeric Backlit LCD Display
- External Adapter, 100-120VAC 150mA, 220-240VAC 100mA, 50/60Hz
- CE, CSA/UL and FCC approved

OPTIONAL ANCILLARIES

- Thermal Printer
- Impact Printer
- Printer Cable
- Ethernet Interface Connector
- Computer Software
- Dust Cover
- Anti Theft Device

| WEIGHTS & DIMENSIONS: DISCOVERY PRO |      |
|-------------------------------------|------|
| Net Weight (kg)                     | 10.2 |
| Width (cm)                          | 20   |
| Depth (cm)                          | 46   |
| Height (cm)                         | 30   |



## TEST SAMPLE PREPARATION

# CNC1 SAMPLE PROFILE CUTTER

Designed and manufactured by Ray-Ran the **Model 1 CNC Sample Profile Cutter** is rapidly becoming the best bench top milling machine in its class.

The 3-axis rapid prototyping vertical milling machine is ideal for cutting hard dense polymer sheets and laminates up to 40mm thick as well as polyethylene and polypropylene pipes that are used within the gas and water industry. The versatility of the machine makes it ideal for Research and Development departments, Laboratories and Universities.

The **CNC1** comes with a bed size of 500mm x 440mm which has tennon slots for clamping down the sample. Travel distances on the X, Y and Z axis are 300mm x 250mm x 100mm respectively and are fitted with high quality lead screws and linear guide rails for precision and accuracy. The high-speed spindle is variable from 2400 rpm to 24000 rpm and has a milling cutter capacity of up to 10mm diameter. Cutting feed ranges up to 1100mm/min can be achieved for a fast-accurate milling process with each axis driven and controlled by advanced stepper motor technology. Pipe work is supplied as standard for air cooling to reduce the cutter temperature when machining components and the machine can also be fitted with dust extraction. The fully enclosed cabinet has an internal light and an electrical safety interlock to ensure operator safety when the cabinet door is in the open position.



For work holding, a standard clamping kit can be supplied or bespoke clamping fixtures can be designed and manufactured by Ray-Ran to meet customer's requirements.

The **CNC1** is supplied with **Flashcut CNC** software which is a dedicated real-time CNC program for use on a PC or Laptop. It interfaces effortlessly with the machine and uses industry standard numerical control G-Code and Letter addresses for programming. **Pre-written CNC programs are available direct from Ray-Ran** for any shape or size of test sample or to any International test standard.

## SAMPLE PROFILE CUTTER (RR/CNC1)

### TECHNICAL SPECIFICATION

- 300mm x 250mm x 100mm XYZ axis (approximate)
- 500 mm x 440 mm table size
- High Quality 0.5KW Router Spindle
- Variable spindle speeds from 2400 rpm to 24000 rpm
- High quality guide rail and lead screw system
- Positional homing switches
- Total enclosure safety cabinet for user protection
- Internal lighting
- Emergency Stop Button
- Computer interface with cabling and Windows™ **Flashcut CNC** PC software
- Low cost CNC programs available to any International test standard
- 220-240 volts 1 ph 50 Hz
- Uses industry standard G & M codes
- Product user manual
- CE declaration certificate
- 1 year return to base warranty

### OPTIONAL ANCILLARIES

- Tooling and clamping kit
- Pipe clamping fixture
- Single flute tungsten carbide milling cutters
- Pre-written CNC programs
- Dust extraction
- I/O spindle expansion board
- Auto tool change system & tool rack
- CAD/CAM software

| WEIGHTS & DIMENSIONS: RR/CNC1 |     |
|-------------------------------|-----|
| Net Weight (kg)               | 110 |
| Width (cm)                    | 98  |
| Depth (cm)                    | 60  |
| Height (cm)                   | 70  |

TECHNICAL SPECIFICATION

- 600mm x 300mm x 175mm XYZ axis (approximate)
- 910 mm x 480 mm table size
- High Quality 1.5KW Router Spindle
- Variable spindle speeds from 2400 rpm to 24000 rpm
- High quality guide rail and lead screw system
- Positional homing switches
- Total enclosure safety cabinet for user protection
- Internal lighting
- Emergency Stop Button
- Computer interface with cabling and Windows™ **Flashcut CNC** PC software
- Low cost CNC programs available to any International test standard
- 220-240 volts 1 ph 50 Hz
- Uses industry standard G & M codes
- Product user manual
- CE declaration certificate
- 1 year return to base warranty

OPTIONAL ANCILLARIES

- Tooling and clamping kit
- Pipe clamping fixture
- Single flute tungsten carbide milling cutters
- Pre-written CNC programs
- Dust extraction
- I/O spindle expansion board
- Auto tool change system & tool rack

WEIGHTS & DIMENSIONS: RR/CNC2

|                 |     |
|-----------------|-----|
| Net Weight (kg) | 150 |
| Width (cm)      | 123 |
| Depth (cm)      | 83  |
| Height (cm)     | 83  |

CNC2 SAMPLE PROFILE CUTTER

The **Model 2 CNC Sample Profile Cutter** designed and manufactured by Ray-Ran is the ultimate bench top 3 axis vertical milling machine used for Test Sample Preparation and 3D rapid prototyping.

It having set new standards for cutting hard dense polymer sheets and laminates up to 100mm thick as well as polyethylene and polypropylene pipes that are used within the gas and water industry. The versatility of the machine makes it ideal for Research and Development departments, Laboratories and Universities.



Totally enclosed for operator safety the **CNC2** comes with a bed size of 910mm x 480mm which has tennon slots for clamping down the sample. Travel distances on the X, Y and Z axis are 600mm x 300mm x 175mm respectively and each slide way is fitted with high quality lead screws and linear guide rails for precision and accuracy. Cutting feed ranges up to 2500mm/min can be achieved for a fast-accurate milling process with each axis driven and controlled by advanced stepper motor technology. Homing switches are fitted as standard to each axis so machine and program co-ordinates can be set easily on a day to day basis ensuring accuracy and repeatability especially if you are using a static clamping fixture for your machining process.

The compact, high torque 1.5KW high speed router spindle which is fitted has a variable speed range from 2400 rpm to 24000 rpm and has a milling cutter capacity of up to 16mm diameter when fitted into industry standard ER25 collets. The Routers spindle speed is controlled via the control panel dial as standard or the optional spindle I/O expansion board, so it becomes fully functional using the **Flashcut CNC** operating software. With this option,

spindle speeds can be set within the G-Code program file and spindle on/off functions selected using specified M-codes. High quality custom made single flute tungsten carbide cutters ensure a superior finish on any machined sample.


To compliment the machine a full clamping kit can be supplied for standard work holding and machining applications. For more specialised work such as machining thick tensile test pieces from industrial polyethylene pipe samples, bespoke work holding fixtures can be designed, manufactured and supplied by Ray-Ran to meet customer's requirements.

Pipe work is supplied as standard for air cooling to reduce the cutter temperature when machining components and the machine can also be fitted with dust extraction. The fully enclosed cabinet has an internal light and an electrical safety interlock to ensure operator safety when the cabinet door is in the open position.

The **CNC2** is supplied with **Flashcut CNC** operating software which is a dedicated real-time CNC program for use on a PC or Laptop. It interfaces effortlessly with the machine and uses industry standard numerical control G-Code and letter addresses for programming. Pre-written CNC programs are available direct from Ray-Ran for any shape or size of test sample or to any International test standard.

If you require a larger machine than the CNC2, Ray-Ran can design, develop and manufacture bespoke alternatives to customer's individual needs. Feel free to contact us to see if Ray-Ran can meet your requirements.

REAL-TIME CNC CONTROL IN A FLASH



**Flashcut CNC** is a dedicated software package which gives high performance low cost CNC results, and is supplied as standard with both models of the Ray-Ran **CNC Sample Profile Cutter**. Effortlessly integrated between machine and PC, Flashcut CNC uses standard numerical G-Code and Letter Addresses for real time CNC control.

The **Flashcut CNC** software has an easy to use, yet powerful operator control panel that puts all of the machine control you need at your fingertips. It has been meticulously engineered to work seamlessly with 64 and

32 bit versions of Windows based operating systems without changing sensitive operating system parameters. It dynamically displays the tool path while sending real time motion commands to the **CNC's** Signal Generator over a high speed USB connection for smooth, dependable, accurate motion on most PC's or Laptops. The system has been designed to accept industry standard G & M codes ensuring full compatibility with leading professional Computer Aided Design and Computer Aided Manufacturing software (CAD/CAM).

For the experienced user of Computer Aided Design software such as AutoCad, simple 2D drawings can be imported into the software as .DXF files using the **Flashcut CNC's** Import feature. Once simple set up parameters are configured the software will generate a simple G-Code program to machine the part. If you have CAM software, programs generated by your Post Processor can also be imported as .FGC files direct to the **Flashcut** software, the program will then generate the cutter path in the viewports on the main control panel.

For the novice, pre-written CNC programs are also available direct from Ray-Ran for any shape or size of test sample to any International test standard, eliminating the need for CAD/ CAM trained operators or skilled CNC programmers.

TECHNICAL SPECIFICATION

- |  |   |   |
|--|---|---|
| <ul style="list-style-type: none"><li>• Easy to use main screen and control panel</li><li>• Program list box</li><li>• G-code program editor</li><li>• Feedrate override selector</li><li>• Step and continuous run modes</li><li>• Digital readout display inc. Program, Machine</li><li>• Relative and Distance to Go co-ordinates</li><li>• Full 5-axis support</li></ul> | <ul style="list-style-type: none"><li>• Tool life management system</li><li>• Very accurate runtime estimate</li><li>• Configuration wizard</li><li>• Advanced contouring</li><li>• Multi-screen viewports</li><li>• DXF import</li><li>• Tool path generator</li><li>• Cycle program creator including, Face Milling, Pocketing,</li><li>• Engraving, Hole Patterns and Custom</li></ul> | <ul style="list-style-type: none"><li>• Cutter compensation feature</li><li>• Main screen display units (english/metric)</li><li>• All G-Code and M-Code variables catered for</li><li>• G-Code teach box</li></ul> |
|--|---|---|

MIN PC REQUIREMENT

- 1 GHz or faster processor
- 1 GB of RAM
- Windows XP, Vista or 7 either 32 or 64 bit
- At least 1 available USB port



## INJECTION MOULDING (RR/TSMA)

# TEST SAMPLE INJECTION MOULDING

The Ray-Ran **Test Sample Injection Moulding Apparatus** has been specifically designed to produce a wide variety of laboratory test samples such as colour plaques, tensile and impact test specimens as well as small components required for mechanical testing procedures. The unique low-cost moulding machine has a compact design for vertical bench mounting making the apparatus ideal for Research and Development Institutes, Universities, laboratories and small workshops.

The apparatus is pneumatically operated with inlet air pressures up to 10 bar (150 psi) which is adequate for moulding most materials including Polypropylene, Polyethylene and PET with injection pressures up to 450 bar (6500 psi)

Digital temperature controls up to 400°C melts the polymer which is then extruded into the mould. Digital tool block heating is included for more efficient polymer flow into the mould with temperatures up to 200°C to prevent freezing.


The apparatus is supplied with a quick release cam operated tool block which will accept low cost tool inserts machined to any shape that will fit within the standard tool block parameters. This makes the apparatus ideal if different

## TEST SAMPLE PREPARATION

samples are required from the same material with minimal wastage and time. Tool Inserts to suit any international standard, e.g. ASTM, ISO, DIN, JIS, etc., can be supplied for tensile dumbbells, colour plaques, impact bars, etc. For samples that cannot be easily moulded within the standard tool block, optional specialised tooling can be supplied to customers' requirements within the parameters of the machine.

The unique interchangeable cylinder liner and die assembly supplied as standard with the apparatus compliments the simplicity of the equipment during and after the moulding process as it enables the operator to change from one colour or material to another in just a few seconds, increasing production and eliminating the need for purging. Simply remove the cylinder liner being used and replace it with a clean one, wait for it to reach temperature, fill the liner with your new material and continue the moulding process. The used liner can then be cleaned at the operator's convenience using a dedicated cleaning solution or ultrasonic cleaning bath. The 47cm<sup>3</sup> capacity of the cylinder liner ensures multiple moulds can be made saving time and reducing valuable material wastage.

To mould a test sample, polymer is loaded into the top of the liner which is



heated to the required temperature to plasticise the polymer. When the injection piston is released, the cylinder liner nozzle automatically locks into the injection gate of the tool block and the piston extrudes the polymer into the moulding plate insert. After the polymer has been allowed to set for a few seconds, the tool block cam is released, and the moulding plate insert withdrawn. The moulded test sample is then removed, the injection gate is cleaned of any residual polymer and the moulding plate insert replaced and clamped in the tool block ready for the next sample to be moulded. For a more convenient flow of polymer into the cylinder liner, an optional polymer feeder hopper can be supplied so the liner is always topped up after each mould is produced ensuring a countless moulding process.


## INJECTION MOULDING (RR/TSMA)

### TECHNICAL SPECIFICATION


- Pneumatically operated
- Maximum air line pressure: 150 psi/10 bar
- Maximum polymer pressure: 6500 psi/450 bar
- Maximum shot size: 47cm<sup>3</sup>
- Maximum sample sizes 175mm (L) x 45 (W) x 12.7 (T)
- Digital Temperature Control
- Automatic nozzle-to-tool locking device.
- Cam-lock Tool Block
- Quick Change Cylinder, Nozzle & Die
- Visual alarms for tool locking
- 110v 60hz and 240v 50Hz
- Product user manual
- Traceable calibration certificate
- CE declaration certificate
- 1 year return to base warranty

### OPTIONAL ANCILLARIES

- Polymer Feeder Hopper
- Low cost moulding inserts available to all international standards
- Quick Change Cylinder, Nozzle & Die
- Air compressor
- Non standard tool block



| WEIGHTS & DIMENSIONS: RR/TSMA |     |
|-------------------------------|-----|
| Net Weight (kg)               | 55  |
| Width (cm)                    | 52  |
| Depth (cm)                    | 36  |
| Height (cm)                   | 111 |



TEST SAMPLE PREPARATION



# HAND OPERATED TEST SAMPLE CUTTING PRESS

Specifically designed to cut test samples such as tensile bars from relatively flexible sheet, the **Ray-Ran Hand Operated Test Sample Cutting Press** will produce samples from plastics, fabrics, boards and paper. In most cases, the thickness of the material that can be cut will not just depend on the strength of the material, but also on the lateral flexibility of the material to allow the cutter to penetrate down into and pass through the material.

The press is supplied with a cutting base of 30cm x 15cm and a polypropylene cutting board to support the material the sample is to be cut from. Four fixing holes in the base plate enable the press to be securely fastened to a rigid bench or worktop to obtain the optimum cutting load. To ensure maximum performance, the cutter must be set up relative to the cutting board so that at full stroke, the cutter will just depress the surface of the cutting board.

The fixing arbour on the press will accept a 25mm diameter cutting tool shank as standard. Other size fixing arbours can be supplied to customer requirements.

The Ray-Ran Hand Operated Cutting Press can generate a cutting force of up to 6KN depending on the strength of the operator.



HAND CUTTING PRESS (RR/HCP)


TECHNICAL SPECIFICATION

- Hand operated
- 30cm x 15cm sample platform
- Cutting force up to 6KN
- Polypropylene cutting board
- Adjustable cutting arbor with locking handle
- Suitable for rubber, plastic, paper, fabric or laminates
- Product user manual
- CE declaration certificate
- 1 year return to base warranty

OPTIONAL ANCILLARIES

- Test sample cutters to any International standard
- Traceable calibration certificate (cutters only)
- Polypropylene cutting board

| WEIGHTS & DIMENSIONS: RR/HCP |    |
|------------------------------|----|
| Net Weight (kg)              | 15 |
| Width (cm)                   | 30 |
| Depth (cm)                   | 25 |
| Height (cm)                  | 50 |



TECHNICAL SPECIFICATION

- Pneumatically operated up to 10 bar (150 psi)
- Automatically operated
- Polycarbonate guard
- Cutting force up to 50KN
- 25cm x 15cm sample platform
- Polypropylene cutting board
- Adjustable cutting arbor with locking handle
- Suitable for rubber, plastic, paper, fabric or laminates
- Product user manual
- CE declaration certificate
- 1 year return to base warranty

OPTIONAL ANCILLARIES

- Test sample cutters to any International standard
- Traceable calibration certificate (cutters only)
- Polypropylene cutting board

WEIGHTS & DIMENSIONS: RR/PCP

|                 |    |
|-----------------|----|
| Net Weight (kg) | 72 |
| Width (cm)      | 30 |
| Depth (cm)      | 30 |
| Height (cm)     | 63 |

PNEUMATICALLY  
OPERATED  
TEST SAMPLE  
CUTTING PRESS

The **Pneumatically Operated Test Sample Cutting Press** is by far the most effortless and safest way to produce cut samples from materials such as plastics, fabrics, boards and paper. The press is robust and compact using a limited amount of space and is ideal for any laboratory or workshop without the need to be bolted to a bench or worktop.

The cutting process is simply done by pulling the protective polycarbonate screen in a downward direction and letting the press do the rest, generating cutting forces up to 50 KN (depending on the input pressure) to easily cutting through sheet materials such as polypropylene and polyethylene up to 8mm thick.


The press is supplied with a cutting base of 25cm x 15cm and a polypropylene cutting board to support the material the sample is to be cut from. A polycarbonate safety guard which will not allow the machine to operate unless it is in the closed position, ensures operator safety. To deliver maximum performance, the cutter must be set up relative to the cutting board so that at full stroke, the cutter will just depress the surface of the cutting board.



The fixing arbour on the press will accept a 25mm dia cutting tool shank as standard. Other size fixing arbours can be supplied to customer requirements.

In most cases, the thickness of the material that can be cut will not just depend on the strength of the material, but also on the lateral flexibility of the material to allow the cutter to penetrate down into and pass through the material.





## TEST SAMPLE PREPARATION


# TEST SAMPLE CUTTERS

Ray-Ran's **Test Sample Cutters** are manufactured from high grade D3 tool steel to give long-lasting cutting-edge sharpness for repeatable test specimens and more accurate test results.

The cutters are heat treated several times during the manufacturing process to remove machining stresses and

are diamond lapped finished for increased cutting-edge wear resistance. The cutters come complete with a spring-loaded ejection system to remove the cut sample from the cutter and protect the cutting edges. A standard fixing shank to suit the Ray-Ran range of presses is fitted to the cutter but non-standard shanks can be manufactured to suit customer's individual requirements.


Cutters can be made to any size or shape required by the user and to any international test standard such as ASTM, ISO, BS, DIN and JIS to name a few.




## TEST SAMPLE CUTTERS (RR/TSC)

### TECHNICAL SPECIFICATION

- Manufactured from high carbon D3 tool steel
- Heat treated for durability
- Diamond lapped finished
- Spring loaded ejection system
- Any shape or size available
- Specimens to suit all national and international standards such as ISO, BS, ASTM, DIN, etc
- Suitable for rubbers, plastics, paper, fabrics and laminates
- Fixing shanks to suit customers own requirements





TEST SAMPLE PREPARATION

AUTOCYCLE

NOTCHING CUTTER

To produce precise stress concentrated notch profiles in Impact and Charpy specimens look no further than the **Ray-Ran Autocycle Notching Cutter**. A world leader by design, this machine will cut accurate notch profiles to international test standards such as ISO 179, ISO 180 and ASTM D256.

At the heart of the machine is a custom-built adaptive Microprocessor Control Unit with LCD screen ensuring feed and cutter speeds are always selected for optimum performance. Adjustments to cutting parameters are made using the integrated membrane keypad and are stored in the memory feature of the microprocessor. Cutter speed ranges from 350 to 2500 rpm are easily controlled, and with minimum feed rates of 0.06mm/rev cutting parameters for even the toughest of materials are easily achieved.

The anti-vibration guide rail system optimises profile accuracy and the multi-sample loading feature of the vice can accommodate up to 16 test specimens 4mm thick, all accurately cut within approximately 20 seconds increasing productivity. The quick auto-return feature ensures the vice is always in the correct place for loading and unloading the test specimens and the micro height adjustment, notch depth comparator and sample setting gauge ensure that the cut notch depth is accurately controlled. A polycarbonate guard



is fitted for operator safety totally enclosing the working area and provides an electrical interlock, so the machine does not operate when the guard is open.

The Dual Tooth Cutting Wheels used on the Autocycle Notching Cutter are designed to reduce cutter load on the sample being notched. They are tungsten carbide tipped, optically ground and diamond lapped to very close tolerances. The standard cutter wheel profile fitted to the apparatus is a 0.25mm radius V notch.

Cutter wheel profiles of 0.1mm and 1.0mm radius as well as 0.8mm and 2.00mm 'U' notch cutter wheels are also available. Non-standard and special notch requirements can also be provided for.

NOTCHING CUTTER (RR/NC)

TECHNICAL SPECIFICATION


- Advanced dedicated microprocessor control
- Easy to read LCD
- Touch membrane keypad
- Memory feature
- Cutter speed range from 350 rpm to 2500 rpm
- Traverse speed range 0.06 mm/rev to 1.0 mm/rev
- Metric and imperial traverse speeds units
- Anti-vibration linear motion slide
- Multi-sample loading
- Cutter wheel 0.25mm radius included
- Notch depth comparator,
- Notch depth gauge block
- Easy clean removable waste tray
- Polycarbonate safety guard
- Electrical Interlock
- Emergency stop button
- Product user manual
- CE declaration certificate
- 1 year return to base warranty
- Notches to International Test Standards, Including: ASTM D256, D6110, ISO 179, 180, 8256, DIN 53453

OPTIONAL ANCILLARIES

- Dual tooth cutter wheel 0.01mm Rad V notch
- Dual tooth cutter wheel 0.25mm Rad V notch
- Dual tooth cutter wheel 1.00mm Rad V notch
- Dual tooth cutter wheel 0.80mm U notch
- Dual tooth cutter wheel 2.00mm U notch

WEIGHTS & DIMENSIONS: RR/NC

|                 |    |
|-----------------|----|
| Net Weight (kg) | 15 |
| Width (cm)      | 54 |
| Depth (cm)      | 33 |
| Height (cm)     | 25 |




## TEST SAMPLE PREPARATION

# MICRO CUTTING MACHINE FOR THIN MATERIALS

Built with absolute simplicity, the Ray-Ran **Micro Cutting Machine** has been specifically designed to cut thin plastic film and other materials into rectangular strips of varied widths using specifically designed cutting heads. It is Ideal for making Tensile and Tear Test strips to various International Test Standards when conventional die cutting is not suitable or recommended. This apparatus has rapidly become a valuable tool from the Ray-Ran product range and is ideal for all types of flexible packaging, film, foils and paper up to 250µm (0.25mm) thick.

Extremely simple to operate, the adjustable spring-loaded head support slides on hardened and ground steel rail guides using linear roller bearings to maintain rigidity and accuracy when cutting. The cutting head attachment is made from high quality tool steel and chrome plated for rust protection and holds 2 custom-made razor-sharp blades accurately spaced to within 0.01mm. The cutter head blades are easily removed when they are worn and need replacing and no further adjustment to the new blades will be required when fitted. With a minimum cut width of 0.25mm, cutting heads can be supplied to cut single or multiple test samples in one single pass increasing your testing capacity and productivity.



The apparatus is supplied with a recessed base which can accept a cutting base pad such as card or rubber sheet which will help protect the cutting blades during the cutting process. A clamping mechanism is also fitted to firmly support the sample material. As standard the cutting stroke length of the apparatus accommodates A4 size samples, but this can be increased to customer's individual requirements.

To operate the Micro Cutting Machine, insert the material under the cutting rails and position it into the clamping mechanism. Slightly slide the cutting head to a position over the material approximately 10 – 20 mm from its clamped edge. Push the cutter down on top of the material and pull the cutter head towards you across the guide rails. Once at the end of its travel, return the cutter head to its start position. Simply release the clamping mechanism, remove the material and trim the sample to the required length.

## MICRO CUTTING MACHINE (RR/MCM)

### TECHNICAL SPECIFICATION

- Simple to operate film cutter
- Cuts film up to 250µm (0.25mm) thick
- Minimum sample width 250µm (0.25mm)
- Hardened and ground guide rail system
- Linear roller bearing guides
- Spring loaded cutter head assembly
- Interchangeable cutter heads
- Multi cut heads available
- Cutter heads made to customers requirements
- Sample clamping rail
- Cutting stroke length A4 as standard.
- Product user manual
- CE declaration certificate
- 1 year return to base warranty

### OPTIONAL ANCILLARIES

- Cutting stroke length to customers requirements
- Cutter heads made to customers requirements
- Multiple blades available on one cutter head
- Cutting Blades (Pack of 10)

| WEIGHTS & DIMENSIONS: RR/MCM |    |
|------------------------------|----|
| Net Weight (kg)              | 5  |
| Width (cm)                   | 40 |
| Depth (cm)                   | 20 |
| Height (cm)                  | 10 |



## SAMPLE & COMPONENT TESTING



## ADVANCED FRICTION TESTER (RR/FT)

# ADVANCED STATIC AND DYNAMIC FRICTION TESTER

The **Advanced Static and Dynamic Friction Tester** brings together accurate load cell measurement with advanced microprocessor technology to give the most accurate, user friendly Friction Tester system to be designed by Ray-Ran. Frictional properties of plastic film and other materials are accurately calculated with specified test conditions selected such as test sled weight, test velocity and temperature making this apparatus ideal for all types of flexible packaging, foils, paper, board, rubber, linoleum, textiles and leather.

Accurate and repeatable test results for Static and Dynamic Friction co-efficient values, with normal separating loads of up to 1 kg, are easily determined for any combination of materials using a standard test sled or specific peel attachment.

To calculate the static friction, the apparatus records the peak friction value at the precise moment that movement occurs after which, the dynamic friction is immediately observed. To measure this peak friction value, the advanced microprocessor continually monitors the output signals of the load cell. Load cell readings are taken beyond the peak load in order to confirm that the slope of the load curve is reducing. When the static frictional value has been determined by the microprocessor the dynamic friction is measured twenty times during the selected travel distance. The result given is the average of these twenty readings. However, by displaying the results graphically on a computer, each of the twenty

readings can be inspected to examine any variation in the test results over the selected travel distance.

The large on-board liquid crystal display (LCD) provides simple on-screen instruction and using the alpha/numeric membrane keypad test parameters such as usernames, material reference numbers and batch numbers are entered and stored in operator lists for future recall and results presentation. Testing parameters such as test velocity, sled weight and settling time are also entered to give the required testing conditions. If Temperature is important to your testing requirement then the Static and Dynamic Friction Tester can be supplied with an optional heated bed. After each test has been conducted the results are displayed on the LCD giving the static and dynamic co-efficient as well as the test load result in grams. Batch statistics are automatically updated after each test and results are shown in graphical and tabular format for analysis when downloaded to



the supplied PC software. For quality control purposes the Static & Dynamic Friction Tester can be set with High and Low limits which are clearly defined when test results are downloaded showing the user instantly if the material is a pass or fail.

The apparatus is supplied as standard with Ray-Ran's **Techni-Test** Windows based PC software for connection to the apparatus via a mini USB or Ethernet connection. Test results can be printed directly from the **Techni-Test software** or can be saved and exported as .CSV files which can be opened with Microsoft Excel in tabular form which can be manipulated to your specific requirements such as preparing a graphical presentation or copying the results to a Word document for a report presentation if required.

Optional vacuum bed is available. Freezer bed is no longer available.

TECHNICAL SPECIFICATION

- Advanced dedicated microprocessor control,
- Touch membrane keypad
- Easy to read liquid crystal display,
- Sequence logic menu auto prompt selection
- Mini USB port
- Ethernet Interface connector for LAN Networking
- Tabular and graphical statistical analysis
- Temperature display
- Variable speed sled velocity up to 1800 mm/minute,
- Test travel distance up to 350mm
- 1 kg load cell
- Sleds to ASTM or ISO test standard supplied
- 180 peel attachment
- Sample cut out template
- Calibration kit included
- Available in either 220-240v 50hz or 110v 60hz
- Fully traceable certificate of calibration,
- Product user manual,
- CE mark of conformity.
- 1 year return to base warranty

TESTING FEATURES

- Operator list
- Material reference list
- Batch ref with data input
- Selectable test type
- Variable sled weight option
- Variable sled velocity option
- Variable test travel distance
- Variable sample settling time
- User defined High and low limits
- Load cell pre-load indicator
- Temperature selector (heated bed option)
- Conforms to: ISO8295, ASTM D1894

OPTIONAL ANCILLARIES


- Heated bed option
- Freezer bed option
- Larger load cells
- Sleds to customer's requirements
- Thermal Printer

| WEIGHTS & DIMENSIONS: RR/FT |    |
|-----------------------------|----|
| Net Weight (kg)             | 24 |
| Width (cm)                  | 67 |
| Depth (cm)                  | 40 |
| Height (cm)                 | 37 |


FRICTION ANALYSIS WITH TECHNI-TEST

**Techni-Test** is an easy to use software package supplied with the **Advanced Static & Dynamic Friction Tester** which allows user defined test data to automatically download to the apparatus for report information and for the operator to analyse all aspects of the friction test.

From the graph, each load cell reading is clearly identified giving accurate frictional data analysis of the materials under test. Multiple tests are clearly visible on the graph for results comparison within the batch and are highlighted in different colours for viewing. Placing the cursor over each plotted point the frictional value can be read on the screen. For internal Quality Control procedures, at a glance high and low limits are clearly displayed on the graph in red on showing instantly if the material is a pass or fail. Tabular results are shown as Test1, Test 2, etc and can be seen by selecting the appropriate tab for a test. Test reports can be printed from the main screen when each test has been completed by highlighting the specific test tab and selecting the **Print Test Report** dialogue box to print the report. Batch statistics such as **Mean** and **Co-efficient of Variation** (COV) are also displayed and are updated after each test result is downloaded.



In Test Results Viewer mode users have the ability to upload saved results from previous test for Material Comparison, Data Manipulation or File Export. Abnormal results that could be caused by static or dust are clearly identified and can be removed from the test data bringing the batch statistics into a normal range ensuring that the test procedure does not have to be repeated saving time and material. Exporting the results file in Viewer Mode is simple. The Export file format is .CSV and can be opened with Microsoft Excel.




# MANUAL BASIC FALLING DART IMPACT TESTER

The Ray-Ran **Manual Basic Falling Dart Impact Tester** provides a basic method for determining the energy that causes thin flexible materials such as plastic film, paper and composite sheets to fail under specified conditions of impact from a free-falling dart. The energy required to cause Initial Failure to the test specimen is expressed in terms of the darts weight which would result in the failure of 50% of specimens being tested. The instrument complies with the ASTM D1709 Method A & B and ISO7765 part 1.

The methods of holding the test specimen and releasing the dart are designed so that they may be accurately repeated for each test performed. Test specimens are simply clamped in the unique two-piece clamp system which gives a constant tension across the surface area of the specimen ensuring test repeatability. A dart of a known weight is then allowed to free fall and hit the sample. The addition or removal of known weights determine the Initial failure of the sample.

To comply with the relevant test standards, a fixed dart drop height of 660mm and 1 ½" spherical radius tup (method A,) and optional 1524mm height and 2" spherical radius (method B) are used. The dart release



mechanism is pneumatically actuated and requires the use of both hands to release the dart to ensure maximum operator safety. A clamp guard is also supplied to ensure the operators safety from falling and rebounding darts which can be removed for easy dart retrieval.

The apparatus comes supplied with digital counters to record Fail and non fails and keep track of the total number of tests conducted.

Fails and Non fails are manually plotted on graph paper to give the test result.

An optional floor mounted stand and method B testing kit are available.

## MANUAL FALLING DART IMPACT TESTER (RR/FDT-MB)

### SPECIFICATIONS (SUPPLIED AS STANDARD)

- Conforms to ISO 7765 and ASTM D1709 methods A
- Bench mounted
- Pneumatically operated
- Simple set up and operation
- Low maintenance.
- Unique constant tension sample clamping
- 2 hand safety operation for clamp and dart release
- Fixed height assembly 0.66mtrs
- Method A 1.5" aluminium tup and weight set.
- Set of 2 x 5g, 8 x 15g, 8 x 30g, 8 x 60g weights
- Fail, non-fail test counters
- Optional ancillaries
- Method B test kit (includes) - 2" stainless steel tup, weight set 2 x 15g, 8 x 45g, 8x 90g, 1524mm height post, remote dart clamp switch.
- Floor mounted stand
- Phenolic, Stainless steel, Aluminium, Composite Darts
- Custom Darts on request

| WEIGHTS & DIMENSIONS: RR/FDT-MB |   |
|---------------------------------|---|
| Net Weight (kg)                 | 50                                      |
| Width (cm)                      | 53                                      |
| Depth (cm)                      | 72                                      |
| Height (cm)                     | 110<br>260 (with optional method B kit) |
| Voltage: 240V 60 Hz / 110v 50Kz |   |

MANUAL ADVANCED FALLING DART  
IMPACT TESTER (RR/FDT-MA)

SPECIFICATIONS  
(SUPPLIED AS STANDARD)

- Conforms to ISO 7765 and ASTM D1709 methods A
- Built in Microprocessor
- Bench mounted
- Pneumatically operated
- Simple set up and operation
- Low maintenance.
- Unique constant tension sample clamping
- 2 hand safety operation for clamp and dart release
- Fixed height assembly 0.66mtrs
- Method A 1.5” aluminium tup and weight set.
- Set of 2 x 5g, 8 x 15g, 8 x 30g, 8 x 60g weights
- Fail, non-fail test counters
- Optional ancillaries
- Method B test kit (includes) - 2” stainless steel tup, weight set 2 x 15g, 8 x 45g, 8x 90g, 1524mm height post, remote dart clamp switch.
- Floor mounted stand
- Phenolic, Stainless steel, Aluminium, Composite Darts
- Custom Darts on request

| WEIGHTS & DIMENSIONS: RR/FDT-MA |   |
|---------------------------------|---|
| Net Weight (kg)                 | 50                                      |
| Width (cm)                      | 53                                      |
| Depth (cm)                      | 72                                      |
| Height (cm)                     | 110<br>260 (with optional method B kit) |
| Voltage: 240V 60 Hz / 110v 50Kz |   |

MANUAL ADVANCED FALLING  
DART IMPACT TESTER

The Ray-Ran Manual Advanced Falling Dart Impact Tester provides simpler method for determining the energy that causes thin flexible materials such as plastic film, paper and composite sheets to fail under specified conditions of impact from a free-falling dart. The energy required to cause Initial Failure to the test specimen is expressed in terms of the darts weight which would result in the failure of 50% of specimens being tested. The instrument complies with the ASTM D1709 Method A& B and ISO7765 part 1.

Utilizing on-board microprocessor technology, the Ray-Ran Microprocessor Controlled Falling Dart Impact Tester has taken falling dart testing machines to an advanced level at an affordable price. Built with operator simplicity in mind, its ease of operation and high accuracy makes it ideal for product development and quality control within production, research & development labs and teaching institutions.

The large on-board liquid crystal display (LCD) provides simple on-screen instruction and using the alpha/numeric membrane keypad test parameters such as usernames, material reference numbers, batch numbers and dart types are entered and stored in operator lists for future recall and results presentation. Selecting the type of test is done at the press of a button and the testing parameters such as drop height; dart weight and sample size are entered to give the required testing conditions.


To comply with the relevant test standards, a fixed dart drop height of 660mm and 1 ½” spherical radius tup (method A,) and optional 1524mm height and 2” spherical radius (method B) are used. The dart release mechanism is pneumatically actuated and requires the use of both hands to release the dart to ensure maximum operator safety. A clamp guard is also supplied to ensure the operators safety from falling and rebounding darts which can be removed for easy dart retrieval.



The methods of holding the test specimen and releasing the dart are designed so that they may be accurately repeated for each test performed. Test specimens are simply clamped in the unique two-piece clamp system which gives a constant tension across the surface area of the specimen ensuring test repeatability. A dart of a known weight is then allowed to free fall and hit the sample. The addition or removal of known weights determine the Initial failure of the sample.

Fails and Non fails are stored and automatically plotted then calculated to give the result making the test easier for the operator. An optional floor mounted stand and method B testing kit are available. The apparatus is supplied as standard with Ray-Ran’s Techni-Test Windows based PC software for connection to the apparatus via a mini USB or ethernet connection. Test results can be printed directly from the Techni-Test software or can be saved and exported as .CSV files which can be opened with Microsoft Excel in tabular form which can be manipulated to your specific requirements such as preparing a graphical presentation or copying the results to a Word document for a report presentation if required.





# ADVANCED FALLING DART IMPACT TESTER A1

Utilizing on-board microprocessor technology, the Ray-Ran Microprocessor Controlled Falling Dart Impact Tester has taken falling dart testing machines to a new level. Built with operator simplicity in mind, its ease of operation and high accuracy makes it ideal for product development and quality control within production, research and development labs and teaching institutions.

The total energy that causes failure of thin flexible materials such as plastic film, paper and composite sheets to fail under specified conditions of impact from a free-falling dart has never been easier to establish. Complying with ASTM D4272, the apparatus uses photo electric cells to measure the Total Energy (joules) to break or cause failure to the sample being tested by measuring the kinetic energy lost by free falling dart.

For perfect energy readings on new samples, a quantity of film will be required to firstly establish the working range of the drop weight required as the dart must completely fall through the sample to give an energy reading during a test. This test is very simple to conduct using the Break Mass test. Once the weight to pass through the sample has been established a

## SAMPLE & COMPONENT TESTING

working range of 50% of the dart weight should be good enough to produce repeatable results. If the break mass is already known, then the Break Energy test will give you your energy reading in Joules

The methods of holding the test specimen and releasing the dart are designed so that they may be accurately repeated for each test performed. Test specimens are simply clamped in the unique two-piece pneumatic clamp system which gives a constant tension across the surface area of the specimen ensuring test repeatability. A dart of a known weight is then allowed to free fall and hit the sample. To comply with the relevant test standard, a fixed dart drop height of 660mm is used and a 1 ½" spherical stainless-steel dart with masses 0.5lb to 2.5lb is used. The dart release mechanism is pneumatically operated for easy release and requires the use of both hands to release the dart to ensure maximum operator safety. A clamp guard is also supplied to ensure the operators safety from falling and rebounding darts and a dart retrieval drawer ensures easy access to the dart after each test.

The large on-board liquid crystal display (LCD) provides simple on-screen instruction and using the alpha/numeric membrane keypad test parameters such as usernames, material reference numbers, batch numbers and dart types are entered and stored in operator lists for future recall and



results presentation. Selecting the type of test is done at the press of a button and the testing parameters such as drop height; dart weight and sample size are entered to give the required testing conditions. Batch statistics of mean, standard deviation and co-efficient of variation are automatically updated after each test and results are shown in graphical and tabular format for analysis when downloaded to the supplied Techni-Test PC software. For quality control purposes High and Low limits can be defined when test results are downloaded showing the user instantly if the material is a pass or fail.

The apparatus is supplied as standard with Ray-Ran's Techni-Test Windows based PC software for connection to the apparatus via a mini USB or ethernet connection.

Test results can be printed directly from the Techni-Test software or can be saved and exported as .CSV files which can be opened with Microsoft Excel in tabular form which can be manipulated to your specific requirements such as preparing a graphical presentation or copying the results to a Word document for a report presentation if required.

The apparatus is supplied with a floor mounted stand only.

### SPECIFICATIONS (SUPPLIED AS STANDARD)

- Conforms to ASTM D4272
- Built in Microprocessor
- Floor Mounted
- Pneumatically operated
- Simple set up and operation
- Low maintenance.
- Unique constant tension sample clamping
- 2 hand safety operation for clamp and dart release
- Fixed height assembly 0.66mtrs
- Method A 1.5" Stainless Steel Tup and weight set.
- Set of 4 x 227g weights

| WEIGHTS & DIMENSIONS: RR/FDT-A1 |     |
|---------------------------------|-----|
| Net Weight (kg)                 | 75  |
| Width (cm)                      | 53  |
| Depth (cm)                      | 72  |
| Height (cm)                     | 260 |
| Voltage: 240V 60 Hz / 110v 50Kz |     |

ADVANCED FALLING DART  
IMPACT TESTER (RR/FDT-A2)

SPECIFICATIONS  
(SUPPLIED AS STANDARD)

- Conforms to ASTM D4272
- Built in Microprocessor
- Floor Mounted
- Pneumatically operated
- Simple set up and operation
- Low maintenance.
- Unique constant tension sample clamping
- 2 hand safety operation for clamp and dart release
- Fixed height assembly 0.66mtrs
- Method A 1.5" Stainless Steel Tup and weight set.
- Set of 4 x 227g weights

WEIGHTS & DIMENSIONS: RR/FDT-A2

|                                 |     |
|---------------------------------|-----|
| Net Weight (kg)                 | 75  |
| Width (cm)                      | 53  |
| Depth (cm)                      | 72  |
| Height (cm)                     | 260 |
| Voltage: 240V 60 Hz / 110v 50Kz |     |




SAMPLE & COMPONENT TESTING

ADVANCED FALLING DART  
IMPACT TESTER A2

Utilizing on-board microprocessor technology, the Ray-Ran Microprocessor Controlled Falling Dart Impact Tester has taken falling dart testing machines to an advanced level at an affordable price. Built with operator simplicity in mind, its ease of operation and high accuracy makes it ideal for product development and quality control within production, research and development labs and teaching institutions.

The break energy that causes thin flexible materials such as plastic film, paper and composite sheets to fail under specified conditions of impact from a free-falling dart has never been easier to establish. Complying with the ASTM D1709 A&B and equivalent ISO standard and ASTM D4272 and quick pass/fail test. The apparatus uses photocells to measure the energy (joules) to break or cause failure to the sample being tested. Usually manual testing requires at least sixty drops for each sample to gain a good result. For the Ray-Ran Microprocessor Controlled Apparatus this is greatly reduced. For perfect energy readings on new samples, a quantity of film will be required to firstly establish the working range of the drop weight required as the dart must completely fall through the sample to give an energy reading during a test. This test is very simple to conduct using the Break Mass test. Once the weight to pass through the sample has been established a working



range of 50% of the dart weight should be good enough to produce repeatable results. If the break mass is already known, then the Break Energy test will give you your energy reading in Joules

The methods of holding the test specimen and releasing the dart are designed so that they may be accurately repeated for each test performed. Test specimens are simply clamped in the unique two-piece pneumatic clamp system which gives a constant tension across the surface area of the specimen ensuring test repeatability. A dart of a known weight is then allowed to free fall and hit the sample. To comply with the relevant test standards, a telescopic shaft for 660mm and 1.5mtr drop heights is used. The dart release mechanism is solenoid actuated for easy release and requires the use of both hands to release the dart to ensure optimum operator safety. A fully electrical interlocked guard is also supplied to ensure the operators safety from falling and rebounding darts and for ruptured samples a dart return shoot is supplied for easy dart retrieval after the test.

SAMPLE & COMPONENT TESTING

The large on-board liquid crystal display (LCD) provides simple on-screen instruction and using the alpha/numeric membrane keypad test parameters such as usernames, material reference numbers, batch numbers and dart types are entered and stored in operator lists for future recall and results presentation. Selecting the type of test is done at the press of a button and the testing parameters such as drop height; dart weight and sample size are entered to give the required testing conditions. At start of each batch of tests a simple calibration of the apparatus is conducted. Batch statistics of mean, standard deviation and co-efficient of variation are automatically updated after each test and results are shown in graphical and tabular format for analysis when downloaded to the supplied Techni-Test PC software. For quality control purposes High and Low limits can be defined when test results are downloaded showing the user instantly if the material is a pass or fail.

The apparatus is supplied as standard with Ray-Ran's Techni-Test Windows based PC software for connection to the apparatus via a mini USB or ethernet connection. Test results can be printed directly from the Techni-Test software or can be saved and exported as .CSV files which can be opened with Microsoft Excel in tabular form which can be manipulated to your specific requirements such as preparing a graphical presentation or copying the results to a Word document for a report presentation if required. Floor mounted only.



ADVANCED FALLING DART  
IMPACT TESTER (RR/FDT-A2)


TECHNI-TEST MAKES A STRONG IMPACT

**Techni-Test** is the easy to use software package supplied as standard with the **Advanced Falling Dart Impact Tester** which allows user defined test data and test results to automatically download from the apparatus for results presentation and data reporting.

Test data downloaded to **Techni-Test** includes operator, material reference, batch reference, test type, impact velocity, hammer weight and impact energy.


From the results graph each impact test value is clearly identified giving accurate data analysis of the material under test. The graphical impact results can be displayed in KJ/M or KJ/M<sup>2</sup> by simply selecting the appropriate results profile tab at the top of the graph. By placing the cursor over each plotted point the impact value can also be read on the screen. For internal Quality Control procedures at a glance high and low limits are clearly displayed on the graph in red showing instantly if the material is a pass or fail. High and low limits are set on the main screen within the **Techni-Test** software so they can be altered easily without having to re-enter test data and re-test material samples. Tabular results are also displayed under the test results tab for quick results processing and comparison analysis between test samples. **Batch** statistics such as **Mean** and **Co-efficient of Variation (COV)** are also displayed and are updated after each test result is downloaded.

Test results can be printed directly from the **Techni-Test** software in a report format or can be saved and exported as .CSV files which can be opened with Microsoft Excel in



tabular form which can be manipulated to your specific requirements such as preparing a graphical presentation or copying the results to a Word document for a report presentation if required

In Test Results Viewer mode, users have the ability to upload saved results from previous test for Material Comparison, Data Manipulation or File Export. Abnormal results that could be caused by air voids are clearly identified and can be removed from the test data bringing the batch statistics into a normal range ensuring that the test procedure does not have to be repeated saving time and material. Test results displayed in viewer mode can also be exported as .CSV files easily.



# ADVANCED UNIVERSAL PENDULUM IMPACT TESTER

The Ray-Ran **Advanced Universal Pendulum Impact Tester** utilises advanced microprocessor technology to determine the energy required to break or rupture specimens such as plastics, composites, ceramics and nonferrous metals to International testing methods for Izod, Charpy and Tension Impact Testing. Along with Pipe Testing, Component Testing and Puncture Impact testing the Ray-Ran Universal Pendulum Impact Tester should more than meet all your testing requirements.

Built with operator simplicity in mind, its ease of operation and high accuracy makes it ideal for product development and quality control within production, research and development labs and teaching institutions. The apparatus is extremely versatile. Impact velocities are variable from 1.5 m/s up to 3.8 m/s and combined with highly accurate variable weight hammers which are electronically released by a solenoid, impact energy ranges up to 50 Joules can be reached.

To calculate the impact energy, the apparatus uses advanced and highly accurate rotary encoder technology to record the lost angle of the hammer after impacting a test sample against the angle recorded after calibrating the hammer. The resulting lost angle is calculated as the test sample's impact energy by the advanced microprocessor system.

The large on-board liquid crystal display (LCD) provides

simple on-screen instruction and using the alpha/numeric membrane keypad test parameters such as usernames, material reference numbers and batch numbers are entered and stored in operator lists for future recall and results presentation.

Selecting the type of test is done at the press of a button and the testing parameters such as impact hammer weight, impact velocity and sample size are entered to give the required testing conditions. At start of each batch of tests, a simple calibration of the apparatus is conducted for bearing resistance and windage and after each test sample has been tested, the results are displayed on the LCD giving the Impact energy in KJ/M, KJ/M<sup>2</sup>, ft/lb.in and ft/lb.in<sup>2</sup>.

Batch statistics of mean, standard deviation and co-efficient of variation are automatically updated after each test and results are shown in graphical and tabular format for analysis when downloaded to the supplied **Techni-Test** PC software. For quality control purposes High and Low limits can be defined when test results are downloaded showing the user instantly if the material is a pass or fail.



- Advanced dedicated microprocessor control,
- Touch membrane Alpha / Numeric keypad
- Easy to read liquid crystal display,
- Sequence logic menu auto prompt selection
- RS232 interface connector
- Ethernet Interface connector for LAN Networking
- Self-calibration procedure for wind and bearing resistance.
- High resolution positional encoder
- Variable pendulum velocity up to 3.8 m/s
- Hammer energies up to 50 Joules
- Results in KJ/M and KJ/M<sup>2</sup>
- Izod, Charpy, Tension, Component, Pipe and Puncture testing.
- Electronic levelling device
- Safety guard
- Solenoid pendulum release with audible pre-warning
- Tabular and graphical statistical analysis
- **Techni-Test** PC software
- Product User Manual
- CE declaration certificate
- Traceable Calibration Certificate
- 1 year return to base warranty

## OPTIONAL ANCILLARIES

- Izod Vice
- Charpy Rest
- Tension Impact Rest (Specimen in bed method)
- Puncture Impact Fixture
- Pipe Rest
- Component fixture to customer requirements
- Variable weight impact hammers up to 50 Joules
- Thermal printer
- Full safety enclosure
- Low temperature chamber


| WEIGHTS & DIMENSIONS: RR/IMT |     |
|------------------------------|-----|
| Net Weight (kg)              | 110 |
| Width (cm)                   | 63  |
| Depth (cm)                   | 50  |
| Height (cm)                  | 76  |

# ADVANCED UNIVERSAL PENDULUM IMPACT TESTER

## TESTING METHODS AVAILABLE



## IZOD TEST

For the Izod test, the specimen is supported and impacted in a cantilevered position with one end of the specimen clamped into the vice. A single swing of the impact hammer is required to break the specimen and the subsequent impact strength calculated. Notched or Un-notched specimens can be tested using the Izod method. Complies to International Test Standards ASTM D256, ASTM D4812 and ISO 180. Impact hammer energies are available from 0.5 J up to 25 J




## CHARPY TEST

The Charpy test requires a standard notched or un-notched test specimen to be supported as a horizontal beam and is impacted on either its narrow face or the width of its cross-section. A single swing of the impact hammer is required to break the specimen and the subsequent impact strength calculated. Complies to International Test Standards ISO 179, ASTM D256, ASTM D6110, DIN 53453, DIN 53753, and BS 7413. Impact hammer energies are available from 0.5 J up to 50 J



## TENSION IMPACT

Tensile impact tests are conducted on material samples that are either too thin or exhibit a high elongation before fracture. The test sample, usually in the shape of a tensile bar, is positioned in a horizontal plane with one end supported to the base of the apparatus and the other end supported by a movable clamp. When the pendulum is released it strikes the movable clamp and breaks the test piece along its horizontal axis, thus applying a pure tensile load to the test sample. Complies to International Test Standard ISO 8256 Method A. Impact hammer energies are available from 0.5 J up to 50 J.



## PIPE TESTING

The test is used for the examination of pipe sections and tubing for impact strength properties as either complete segments or small sections of pipes. Sample diameter dimensions up to 25 mm can be tested with hammer energies of 7.5 - 15 J. The test specimen is supported as a horizontal beam similar to Charpy tests conditions and complies to



## SAMPLE & COMPONENT TESTING

International test Standards ISO 7628 and ISO 9854. Sample diameter dimensions up to 25 mm (0.98 in) can be tested with hammer energies of 7.5 - 15 J (5.6 - 11.1 ft-lbs) or 50 J (36.9 ft-lbs), as defined in the ISO standards.

## COMPONENT TEST

The component test examines the influences of both the design and the manufacturing processes of a component by measuring the impact energy in Joules, required to break a complete or part of a section of a finished component. If a component has a critical cross-sectional area that needs to be tested, then by entering this area into the microprocessor the energy in KJ/M<sup>2</sup> can be obtained.




The apparatus is supplied as standard with Ray-Ran's **Techni-Test** Windows based PC software for connection to the apparatus via a mini USB or Ethernet connection. Test results can be printed directly from the **Techni-Test** software or can be saved and exported as .CSV files which can be opened with Microsoft Excel in tabular form which can be manipulated to your specific requirements such as preparing a graphical presentation or copying the results to a Word document for a report presentation if required. For quality control purposes, High and Low limits can be defined when test results are downloaded showing the user instantly if the material is a pass or fail

## TECHNI-TEST MAKES A STRONG IMPACT

**Techni-Test** is the easy to use software package supplied as standard with the **Advanced Universal Pendulum Impact Tester**, which allows user defined test data and test results to automatically download from the apparatus for results presentation and data reporting.

Test data downloaded to **Techni-Test** includes operator, material reference, batch reference, test type, impact velocity, hammer weight and impact energy.



(COV) are also displayed and are updated after each test result is downloaded.

From the results graph each impact test value is clearly identified giving accurate data analysis of the material under test. The graphical impact results can be displayed in KJ/M or KJ/M<sup>2</sup> by simply selecting the appropriate results profile tab at the top of the graph. By placing the cursor over each plotted point the impact value can also be read on the screen. For internal Quality Control procedures at a glance high and low limits are clearly displayed on the graph in red showing instantly if the material is a pass or fail. High and low limits are set on the main screen within the **Techni-Test** software so they can be altered easily without having to re-enter test data and re-test material samples. Tabular results are also displayed under the test results tab for quick results processing and comparison analysis between test samples. **Batch** statistics such as **Mean** and **Co-efficient of Variation**

Test results can be printed directly from the **Techni-Test** software in a report format or can be saved and exported as .CSV files which can be opened with Microsoft Excel in tabular form which can be manipulated to your specific requirements such as preparing a graphical presentation or copying the results to a Word document for a report presentation if required

In Test Results Viewer mode users have the ability to upload saved results from previous test for Material Comparison, Data Manipulation or File Export. Abnormal results that could be caused by air voids are clearly identified and can be removed from the test data bringing the batch statistics into a normal range ensuring that the test procedure does not have to be repeated saving time and material. Test results displayed in viewer mode can also be exported as .CSV files easily.

ADVANCED UNIVERSAL PENDULUM  
IMPACT TESTER WITH PUNCTURE IMPACT

TECHNICAL SPECIFICATION

- Advanced dedicated microprocessor control,
- Touch membrane Alpha / Numeric keypad
- Easy to read liquid crystal display,
- Sequence logic menu auto prompt selection
- RS232 interface connector
- Ethernet Interface connector for LAN Networking
- Self-calibration procedure for wind and bearing resistance.
- High resolution positional encoder
- Variable pendulum velocity up to 3.8 m/s
- Hammer energies up to 50 Joules
- Results in KJ/M and KJ/M<sup>2</sup>
- Izod, Charpy, Tension, Component, Pipe and Puncture testing.
- Electronic levelling device
- Safety guard
- Solenoid pendulum release with audible pre-warning
- Tabular and graphical statistical analysis
- **Techni-Test** PC software
- Product User Manual
- CE declaration certificate
- Traceable Calibration Certificate
- 1 year return to base warranty

OPTIONAL ANCILLARIES

- Izod Vice
- Charpy Rest
- Tension Impact Rest (Specimen in bed method)
- Puncture Impact Fixture
- Pipe Rest
- Component fixture to customer requirements
- Variable weight impact hammers up to 50 Joules
- Thermal printer
- Full safety enclosure
- Low Temperature Chamber

WEIGHTS & DIMENSIONS: RR/IMT

|                 |     |
|-----------------|-----|
| Net Weight (kg) | 110 |
| Width (cm)      | 63  |
| Depth (cm)      | 50  |
| Height (cm)     | 76  |



SAMPLE & COMPONENT TESTING

ADVANCED UNIVERSAL PENDULUM IMPACT TESTER  
WITH PUNCTURE IMPACT

Determining the most cost-effective packaging material specification for an application can be a daunting proposition. The impact strength can be crucial if the handling of the package cannot be controlled, especially if objects with sharp profiles are to be packaged and reliable information is not available on the material behaviour under different load conditions. The Ray-Ran **Advanced Universal Pendulum Impact Tester** has been designed to overcome these problems using the optional **Puncture Impact** fixture which makes it ideal for research and development within the packaging industry.

The Puncture Impact Tester can examine the material behaviour of thin brittle or ductile materials in a combination of different ways. By using different Tup Strikers, the mode of failure can be analysed to represent direct comparisons with other materials and produce in-house quality control standards, or to simulate actual service conditions. For instance, if impact elongation properties are required for ductile materials, then a large spherical radius Tup Striker is used to induce a high degree of lateral elongation of the material before the Tup Striker can penetrate through the material to the point of failure. If, however, a more brittle failure (minimum elongation of the material) is required, a more pointed Tup Striker is attached to the hammer to impact on the sample. If required, the Tup Striker can simulate the actual

object to be packaged to see how the packaging material would stand up to poor handling in service.


It is possible that the direction of material flow of the specimen



SAMPLE & COMPONENT TESTING

could influence the test result therefore the test sample takes the form of a circular specimen, up to 500 microns thick which is clamped in a sample holder. To ensure repeatability, a rubber “O” ring is used to exert the same tension across the test sample and to prevent the sample moving when under test. The unique design of the variable weight pendulum hammer enables sheet material samples to be tested and punctured for ductile and brittle failure analysis. An impact hammer energy range of up to up to 25 Joules and testing velocities up to 3.8 m/s can be achieved ensuring the toughest of materials can be tested. To calculate the impact energy, the apparatus records the lost angle of the hammer after impacting a test sample against the angle recorded after calibrating the hammer. The resulting lost angle is calculated as the test sample’s impact energy by the advanced microprocessor system.

The large on-board liquid crystal display (LCD) provides simple on-screen instruction and using the alpha/numeric membrane keypad test parameters such as usernames, material reference numbers and batch numbers are entered and stored in operator lists for future recall and results presentation. Selecting the type of test is done at the press of a button and the testing parameters such as impact hammer weight, impact velocity and sample size are entered to give the required testing conditions. At start of each batch of tests, a simple calibration of the apparatus is conducted for bearing resistance and windage and after each test sample has been tested, the results are displayed on the LCD giving the impact energy in Joules. Batch statistics of mean, standard deviation and co-efficient of variation are automatically updated after each test and results are shown in graphical and tabular format for analysis when downloaded to the supplied **Techni-Test** PC software. For quality control purposes High and Low limits can be defined when test results are downloaded showing the user instantly if the material is a pass or fail.





FALLING WEIGHT IMPACT TESTER (RR/FWT)

# UNIVERSAL FALLING WEIGHT IMPACT TESTER

Designed and built with operator simplicity in mind, the Ray-Ran **Universal Falling Weight Impact Tester** is simply the best in its class. It is used to determine the energy required to break or rupture test specimens such as pipe, sheet, laminates, composites, ceramics and nonferrous metals for material and component evaluation to international testing methods such as ISO6603, ISO3127, ASTMD2444 and UL 651 to name a few.

Utilising both pneumatic and electrical functions, the apparatus is extremely versatile. The standard machine is supplied with a 2-meter variable drop height system which can reach impact velocities up to 6.26 m/s. With an impact energy range of up to 314 joules, even the toughest of materials can easily be tested.

Built with operator safety in mind, the apparatus has a full electrical interlock system preventing operation of the apparatus if a safety guard is open. It is supplied as standard with a solenoid operated carrier release mechanism for simple operation and a unique pneumatic carrier return system to ensure the carrier is safely returned to its set drop height for the next test to be performed.

SAMPLE & COMPONENT TESTING



A twin guide rail system is used to ensure a smooth repetitive drop of the load carrier perpendicular to the sample being tested. Load carriers are available for different energy ranges and are fitted with low friction bearing guides which prevent velocity and potential energy being lost as the carrier descends. Impact tups are supplied to international test standard methods, or to customers own individual requirements, which can simulate the mode of failure in actual service conditions or to analyse specific in-house quality control standards.

To conduct a test, prepared samples are placed in the enclosed chamber on their supporting fixture which can be adjusted vertically to accommodate different test sample sizes. On the standard machine pipe diameters of up to 400mm Ø can be tested easily. If the optional pneumatic clamping fixtures are used, the integrated safety system of the apparatus ensures that the clamp does not operate until the cabinet door is in the closed position.

## THE METHODS OF TESTING ON THE UNIVERSAL FALLING WEIGHT IMPACT TESTER ARE

### ROUND THE CLOCK METHOD

Lengths of pipe from a batch or production run are subjected to blows from a known mass which falls from a specified height as per test standard ISO 3127 & EN744.

### STAIRCASE METHOD

Lengths of pipe are subjected to blows from a known mass and shape but is dropped from differing heights depending on the results of each blow as per test standard EN1411.

The **Universal falling Weight Impact Tester** is ideal for product development and quality control within production, research and development labs and teaching institutions and will more than meet your testing requirements. Various drop heights are available as well as larger bases to accept larger pipe diameters and Tup (Striker) shapes to customers individual requirements.

FALLING WEIGHT IMPACT TESTER (RR/FWT)

### TECHNICAL SPECIFICATION

- 2 meter drop height test machine
- Twin guide rail system
- Low friction carrier bearing system
- Solenoid release mechanism
- Pneumatic carrier return
- Mechanical raise/lower fixture platform
- Pipe, sheet, Izod, Charpy, Tension, Component testing compatible.
- Variable velocity up to 6.26 m/s
- Impact energies up to 314 Joules
- Electrical safety interlock
- 110 or 240 volt available
- Operating pressure 100 PSI (7 bar)
- Conforms to ISO6603, ISO3127 and ASTM D2444, etc
- Product user manual
- CE declaration certificate
- Traceable calibration certificate
- 1 year return to base warranty

### OPTIONAL ANCILLARIES

- Extended drop height
- Increased chamber capacity
- "V" rest pipe fixture 12mm – 75mm dia
- "V" rest pipe fixture 75mm – 400mm dia
- Plate or Plaque specimen holder (manual clamp)
- Plate or Plaque specimen holder (Pneumatic clamp)
- Izod Fixture
- Charpy Fixture
- Customized Fixtures
- Anti-rebound device
- Drop Height Power Feed Adjustment
- Digital Height Scale
- Velocity indicator
- Drop weight carriers from 0.5kg – 16kg as standard
- Customized drop weight carrier
- Impact Heads (Tups) to International test standards
- Customized Impact Heads (Tups)
- Variable weight system
- Low Temperature Chamber

| WEIGHTS & DIMENSIONS: RR/FWT          |       |
|---------------------------------------|-------|
| Net Weight (kg)                       | 270*  |
| Width (cm)                            | 950*  |
| Depth (cm)                            | 450*  |
| Height (cm)                           | 2500* |
| *standard 2 meter drop height machine |       |

## AVAILABLE OPTIONS

### ANTI REBOUND

The apparatus can be supplied with an electrically operated anti-rebound system to prevent the load carrier rebounding on to the sample after its initial impact preventing further damage to the test specimen.

### DIGITAL HEIGHT SCALE

This optional feature can be fitted to accurately determine the drop height of your weight carrier. The fitted LED digital readout clearly displays the drop height in mm or inches.

### DROP HEIGHT POWER FEED ADJUSTMENT

To further enhance the apparatus, an optional powered lift can be fitted to raise and lower the drop height settings. This feature is especially useful if your testing drop heights are in excess of 2 meters.

### VELOCITY INDICATOR

Fitted as an option, this device accurately determines the falling velocity of the weight carrier making impact energy calculations even more precise.

# UNIVERSAL FALLING WEIGHT IMPACT TESTER

## FIXTURES


By changing the supports and clamping mechanisms, along with the shape of the Tup, and variable weights the Ray-Ran **Universal Falling Weight Impact Tester** can be manufactured to meet any test specification or costumers requirement making the equipment truly bespoke. Typical test fixtures available include

### PIPE TESTING "V" REST

The test is used for the examination of pipe sections and tubing for impact strength properties as complete segments of pipe. Sample diameters from as little as 12mm Ø up to 400mm Ø can be tested on the standard machine. The pipe is supported in a 120° angled "V" support as is impacted by the load carrier in a single drop. The support rests can be manufactured to any test standard or customer requirement.

### FLAT PLATE / PLAQUE TESTING

Testing of flat plates or plaques such as laminates and



composites can easily be tested using the flat plate clamping fixture. It is offered with various support rings to test to various international test standards and customers own specific testing requirements and can be either manually or pneumatically operated.

Fixtures for conducting **Izod** and **Charpy** tests are also available.

For specific component testing Ray-Ran can design and manufacture bespoke clamping fixtures to meet customers individual testing requirements and parameters.





# ADVANCED HDT/VICAT SOFTENING POINT APPARATUS 2 & 4 STATION

Designed and built to cover multiple international testing standards, the Ray-Ran **Advanced HDT/Vicat Apparatus** utilises microprocessor technology to accurately determine the deflection and softening point characteristics of all thermoplastic test specimens.

The HDT/Vicat apparatus allows for multiple simultaneous testing of samples depending on the amount of test stations available, either 2 or 4. The simple manually operated raise and lower function of the test stations ensure easy access to the test sample supports for sample loading and retrieval before and after each test. Each test station is fitted with a PT100 platinum resistance thermometer which accurately records the temperature next to the sample under test to 0.1°C and an electronic displacement transducer which measures the sample displacement to 0.01mm as standard or to 0.001mm as an option. The on-board microprocessor ensures test result accuracy and repeatability and the built-in liquid crystal display (LCD) provides simple on-screen instructions reducing user error.

Test parameters are easily selected using the onboard membrane keypad. A simple data selection process and yes/no prompts make the operating procedure very simple to undertake. The microprocessor's temperature control function ensures the ramping rates of either 50°C/h or 120°C/h are kept within the specified test standard requirements as well as performing non standard ramp rates to customers own specific needs. Test temperatures of 300°C are easily achieved and to ensure optimum safety at higher temperatures, the Nitrogen Blanket option is recommended.



## THE METHODS OF TESTING ON THE APPARATUS ARE HDT – HEAT DEFLECTION / DISTORTION TEST

A standard sized test specimen is subjected to a bending stress, whilst the temperature is raised at a uniform rate. The temperature at which the specified deflection occurs is measured and recorded. Testing is carried out in accordance with the ISO 75 (parts 1, 2 and 3) and ASTM D648 Test Standards. For this test, a required fibre stress of 0.45, 1.8, or 8.00 MPa is easily selected. Custom fibre stresses can also be managed by the microprocessor if required. The unique binary weight system is used to apply the required fibre stress to the test sample and is automatically calculated by the microprocessor based on the sample size and span supports. Temperature ramp rates, sample size, span and deflection values are also easily entered into the testing parameters of each station. The HDT span supports of 100mm or 64mm are easily adjusted on each tool station to suite your testing method of **Flat wise** or **Edgewise** sample testing. Each machine is supplied with HDT test nibs with 3mm radius which are easily attached to the load displacement rods for testing in accordance with relevant international test standards.

ADVANCED HDT/VICAT SOFTENING  
POINT APPARATUS (RR/HDV2 – RR/HDV4)

TECHNICAL SPECIFICATION

- HDT/Vicat testing enabled
- Manual Raise/Lower of test stations
- Advanced microprocessor control
- 2 or 4 sample test stations
- Digital temperature control
- Temperature range to 300°C
- Oil bath stirrer
- Solenoid operated cooling system
- Integrated safety thermostat
- HDT Heads (1 per station)
- Vicat Nibs (1 per station)
- Standard Fibre Stress 0.45, 1.8 or 8.00 MPa
- Used defined fibre stress for HDT testing available
- Standard Vicat penetration of 0.1mm or 1.00mm
- User defined penetration depth
- for Vicat testing available
- Temp resolution +/- 0.1°C
- Temperature ramp rates of 50° or 120° C/hr
- User defined ramp rates available
- HDT span supports 64mm or 100mm
- Resolution +/- 0.01mm
- Bath capacity 8 litres
- RS 232 connection/RS232 Cable
- Windows File Capture Software
- Fully traceable certificate of calibration
- Product user manual
- CE declaration certificate
- 1 year return to base warranty
- 110v – 240v 16amp
- Conforms to ISO 75-1, 75-2, 75-3, 306 and ASTM D648, D1525, D5944, D5945

OPTIONAL ANCILLARIES

- Displacement Transducer 0.001mm
- Binary HDT Test Weights (1 set required per station)
- 1.00Kg Vicat Test Weight (1 required per station)
- 5.00Kg Vicat Test Weight (1 required per station)
- Light Weight Load rods to test to .45 MPa
- Heat Transfer Medium 5ltrs
- Nitrogen Blanket
- Air to Nitrogen Extractor
- Thermal Printer
- Thermal Printer Paper
- Water Chiller Unit

| WEIGHTS & DIMENSIONS:<br>HDT/VICAT SOFTENING POINT APPARATUS |      |      |
|--|------|------|
|  | HDV2 | HDV4 |
| Net Weight (kg)  | 25   | 40   |
| Width (cm)   | 36   | 46   |
| Depth (cm)   | 53   | 53   |
| Height (cm)  | 50   | 50   |



SAMPLE & COMPONENT TESTING

VICAT (VST) –  
SOFTENING POINT TEST

A specified needle (indenter) penetrates a specified distance into a sample with a specified load, whilst the temperature is raised at a uniform rate. The temperature at which the sample was penetrated is recorded as the Vicat Softening Temperature (VST). Testing is carried out in accordance with the ISO 306 and ASTM D1525 Test Standards.

For this test, penetration loads are easily selected via the microprocessor for use with the single loading weights of 10 N and 50N. The supplied cylindrical


indenter test nib with 1mm<sup>2</sup> surface area is attached to the load displacement rod for sample testing in accordance with relevant testing standard.

The apparatus is supplied as standard with an integrated solenoid operated cooling system which is automatically activated after each test. It can be connected to either a standard water supply or to the optional chiller unit for improved heat reduction. The integrated cooling coil ensures rapid heat loss back to start temperature conditions within a short time period increasing production.

The self-calibration feature of the apparatus ensures that the test results remain accurate. Unique calibration plugs are supplied with each machine which simply plug into the PRT socket for you to carry out the calibration procedure.

Test results are downloaded to the supplied file capture software via the RS232 serial connector and are displayed in tabular format. Results can also be saved as .CSV files which can be opened with Microsoft Excel for data manipulation and report presentation. If required, the optional thermal printer can be supplied and fitted to the apparatus for direct results printout.

Built with operator simplicity in mind, its ease of operation and high accuracy makes the Ray-Ran Advanced HDT/Vicat Apparatus ideal for product development and quality control within production, research and development labs and teaching institutions and should more than meet all your testing requirements.



TECHNICAL SPECIFICATION

- HDT/Vicat testing enabled
- Advanced microprocessor control
- Power rise and fall of test stations
- Dual pump stirrer system
- 6 sample test stations
- Digital temperature control
- Temperature range to 300°C
- Solenoid operated cooling system
- Integrated safety thermostat
- HDT Heads (1 per station)
- Vicat Nibs (1 per station)
- Standard Fibre Stress 0.45, 1.8 or 8.00 MPa
- Used defined fibre stress for HDT testing available
- Standard Vicat penetration of 0.1mm or 1.00mm
- User defined penetration depth for Vicat testing available
- Temp resolution +/- 0.1°C
- Temperature ramp rates of 50° or 120° C/hr
- User defined ramp rates available
- HDT span supports 64mm or 100mm
- Resolution +/- 0.01mm
- Bath capacity 18 litres
- Built in printer
- RS 232 connection & cable
- Windows data logging software
- Fully traceable certificate of calibration
- Product user manual
- CE declaration certificate
- 1 year return to base warranty
- 240v 16amp
- Conforms to ISO 75-1, 75-2, 75-3, 306, and ASTM D648, D1525, D5944, D5945

OPTIONAL ANCILLARIES

- Displacement Transducer 0.001mm
- Binary HDT Test Weights (1 set required per station)
- 1.00Kg Vicat Test Weight (1 required per station)
- 5.00Kg Vicat Test Weight (1 required per station)
- Light Weight Load rods to test to .45 MPa
- Heat Transfer Medium 5ltrs
- Nitrogen Blanket
- Air to Nitrogen Extractor
- Water Chiller unit

WEIGHTS & DIMENSIONS: HDT/VICAT  
SOFTENING POINT APPARATUS (HDV6)

|                 |    |
|-----------------|----|
| Net Weight (kg) | 50 |
| Width (cm)      | 90 |
| Depth (cm)      | 60 |
| Height (cm)     | 40 |

# ADVANCED 6 STATION HDT/VICAT SOFTENING POINT APPARATUS

For maximum testing capability look no further than the Ray-Ran **Advanced 6 Station HDT/Vicat Apparatus**. Built for multiple simultaneous testing of 6 samples to HDT or Vicat testing methods, the apparatus uses dedicated microprocessor technology to accurately determine the deflection and softening point characteristics of all thermo plastic test specimens. Its operator simplicity and high accuracy, makes the apparatus ideal for product development and quality control and has been designed to meet multiple International Testing Standards.

On start up the advanced on board microprocessor starts work immediately by conducting a self-diagnostic procedure to ensure test result accuracy and reliability. The built-in liquid crystal display (LCD) provides simple on-screen instructions and test parameters are easily selected using the onboard membrane keypad. A simple data selection process using yes/no prompts make the setup and operating procedure very simple to undertake and reduces user error and operator training. The microprocessor's temperature control function

ensures the ramping rates of either 50°C/h or 120°C/h are kept within the specified test standard requirements as well as performing nonstandard ramp rates to customers own specific needs

Each test station is fitted with a PT100 platinum resistance thermometer, which accurately records the temperature next to the sample under test to 0.1°C and an electronic displacement transducer, which measures the sample



displacement to 0.01mm as standard or to 0.001mm as an option.

The integrated power raise and lower function of the test stations ensures safe easy access to the test sample supports for sample loading and retrieval before and after each test. For maximum temperature stability across each of the test stations the apparatus is fitted with a dual pump oil stirrer system. Test temperatures of 300°C are easily achieved and to ensure optimum safety at higher temperatures, the Nitrogen Blanket option is recommended.

The apparatus is supplied, as standard with an integrated solenoid operated cooling system which is automatically activated after each test. It can be connected to either a standard water supply or to the optional chiller unit for improved heat reduction. The integrated cooling coil ensures rapid heat loss back to start temperature conditions within a short time period increasing production.

The one touch self-calibration feature of the apparatus ensures that the machine and test results remain accurate. A unique calibration interface unit is supplied with each machine that is simply connected to the apparatus for you to carry out the calibration procedure.

The apparatus is supplied, as standard with advanced data logging software which connects via the RS232 serial interface connector. The software records the temperature/deflection profile of each station in real time during a test. Test results are displayed in tabular format and reports can easily be generated by the software for results presentation. If required, test results can also be saved as .CSV files which can be opened with Microsoft Excel for data manipulation, further report presentation or for importing into customer specific software.





## SAMPLE & COMPONENT TESTING

# COLD FLEX TESTER

Designed and built by Ray-Ran, the **Cold Flex Tester** determines how low temperatures affect the torsional stiffness properties of flexible materials, including polyvinyl chloride extrusion compounds, by measuring the temperature at which a test specimen is twisted through a known angular displacement using a pulley system by a specified torque. The apparatus complies with BS2782: Part 1: Method 104B, ISO458/1, 458/2 and ASTM D1043 International Test Standards.

The apparatus consists of a low temperature bath with electronic PID temperature controller, heater and a PT100 platinum resistance thermometer accurate to 0.1°C which accurately controls the heating cycle of the test. To ensure temperature stability within the bath during the cooling and heating cycle, a stirrer motor system is also fitted to the apparatus. Suitable means should be used to ensure the liquid medium can be lowered to the test temperatures required such as adding solid carbon dioxide (dry ice).

The specimen clamping mechanism is manufactured from stainless steel and has one fixed side with the other side mounted in low friction bearings which can rotate. The torque is applied to the test sample by placing known weights to the pulley system which causes the clamped test specimen to rotate or twist about its vertical axis. Angular deflection measurements (in degrees) are then taken from the integrated radial dial indicator.

The actual test is very simple to conduct. The temperature



of the bath is lowered to a temperature that will cause the sample to approximately twist a known angle of deflection. Heat or increased torque is then applied to the sample to increase its twist to a given angle of deflection subject to the international Test Standard you are working to. From the temperature and the angle of deflection readings taken, the apparent torsional modulus of elasticity can be determined.

Simple and unique, this apparatus has become a valuable addition to the Ray-Ran product range of materials testing equipment.

## COLD FLEX TESTER (RR/CFT)

### TECHNICAL SPECIFICATION


- Measures apparent torsional modulus of elasticity
- Integrated low temperature bath
- PID electronic temperature control
- Resolution 0.1°C
- PT100 PRT sensor accurate to 0.1°C
- Integrated stirrer motor
- Stainless steel clamping mechanism
- Low friction bearing assembly
- Integrated radial dial indicator (degrees)
- Torque load weights supplied
- Product user manual
- CE declaration certificate
- 1 year return to base warranty
- 110v 60Hz – 240v 50Hz
- 10 amp
- Complies with BS2782: Part 1: Method 104B, ISO458/1 & 458/2
- ASTM D1043

### OPTIONAL ANCILLARIES

- Heat transfer medium (silicone oil) 10 ltrs

### WEIGHTS & DIMENSIONS: RR/CFT

|                 |    |
|-----------------|----|
| Net Weight (kg) | 15 |
| Width (cm)      | 30 |
| Depth (cm)      | 45 |
| Height (cm)     | 40 |



## SAMPLE & COMPONENT TESTING

# LOW TEMPERATURE BRITTLENESS APPARATUS

The Ray-Ran **Low Temperature Brittleness Tester** has been designed to study the effects low temperatures have on plastic materials when subject to impact loads and conforms to ASTM D746 and ISO R974 International Test Standards. Many materials which are flexible at normal ambient temperatures become brittle at low temperatures which severely affect the materials impact strength. By plotting the percentage failure rate of the test samples against the reduction in temperature, the brittleness temperature of the test samples can be determined.

The apparatus consists of a low temperature bath with electronic PID temperature controller and a PT100 platinum resistance thermometer accurate to 0.1°C which accurately controls the heating cycle of the test and records the temperature next to the test samples. Precise temperature control is readily achieved by means of a built-in heater which, in conjunction with controlled addition of refrigerant, enables any temperature down to -70 °C to be readily attained and



held. To ensure temperature stability within the bath during the cooling and heating cycle a stirrer motor system is also fitted to the apparatus. Suitable means should be used to ensure the liquid medium can be lowered to the test temperatures required, such as adding solid carbon dioxide (dry ice) to Methyl Alcohol. The tester is fitted with a safety interlocked guard for operator safety.

The test procedure is very simple to conduct, and a sufficient number of test specimens should be tested so the Brittleness Temperature can be calculated on a statistical basis. Test specimens are clamped in a stainless-steel cantilever beam fixture which is immersed into the low temperature bath. The specimens are allowed to condition at a maintained temperature in the heat transfer medium and are then simultaneously subjected to a single blow given by a free-falling pendulum impact hammer through the means of an intermediate striker at a constant velocity of 2.1 m/s. The samples are then removed and examined to determine whether failure has occurred. Temperatures are then increased or decreased to find the temperature at which 50% of the specimens fail. This is defined as the Brittleness Temperature.

## LOW TEMPERATURE BRITTLENESS APPARATUS (RR/LTB)

### TECHNICAL SPECIFICATION

- Low temperature brittleness determination
- Integrated low temperature bath
- PID electronic temp control
- Resolution 0.1°C
- PT100 PRT sensor accurate to 0.1°C
- Integrated stirrer motor
- Temperature range - 70°C to 100°C
- Stainless steel cantilever clamp
- Easy specimen removal
- Manual pendulum release
- Low friction bearing assembly
- Integrated specimen striker head
- Impact velocity 2.1 m/s
- Test sample cutting die
- Product user manual
- CE declaration certificate
- 1 year return to base warranty
- 110v 60Hz – 240v 50Hz
- 10 amp
- Product user manual
- CE declaration certificate
- 1 year return to base warranty
- Complies with: ISO 974 AND ASTM D746

### OPTIONAL ANCILLARIES

- Specimen cutter Type A
- Specimen cutter Type B
- Specimen cutter Type C
- Specimen Clamps to ISO or ASTM test standard.

### WEIGHTS & DIMENSIONS: RR/LTB

|                 |    |
|-----------------|----|
| Net Weight (kg) | 35 |
| Width (cm)      | 60 |
| Depth (cm)      | 25 |
| Height (cm)     | 50 |



## SAMPLE & COMPONENT TESTING

# ENVIRONMENTAL STRESS CRACKING APPARATUS

### The Ray-Ran **Environmental Stress Cracking Apparatus**

determines the susceptibility of ethylene plastics to stress-cracking when exposed to different environments such as soaps, wetting agents, oils and detergents. When the material is under certain conditions of stress, these reagents can accelerate the cracking process which is one of the most common causes of unexpected brittle failure of polymers known. Designed in accordance with the ASTM D1693 international test standard, the apparatus is extremely cost effective and very simple to use.

Manufactured to a high standard, the apparatus is supplied with electronic digital temperature control, temperature bath and all tooling accessories. The integrated temperature bath fitted with digital temperature controller, heater and a PT100 platinum resistance thermometer accurately maintains the test temperature, usually 50 °C or 100 °C ± 0.5 °C. To ensure temperature stability within the bath, a stirrer motor system is fitted, and the built-in electronic timer monitors the duration of the test. The bath is fitted with a stainless-steel sample rack which holds up to 48 glass test tubes which are immersed in



a heat transfer medium such as silicone oil. Rubber stoppers to seal the tubes and test sample specimen holders in either stainless steel or brass (optional) are also supplied for each tube. The specimen holder is designed to maintain a constant stress on the mid-section of the test samples with up to 10 specimens fitted into a specimen holder at one time.

### THE UNIQUE SET OF TOOLING ACCESSORIES REQUIRED TO PRODUCE YOUR TEST SAMPLES INCLUDE:

- **TEST SAMPLE CUTTER** – 38 mm x 13 mm.
- **SAMPLE NICKING JIG** – Used to make a controlled imperfection in the test sample.
- **SAMPLE BENDING CLAMP** – Used to induce a concentrated stress in the sample.
- **SAMPLE TRANSFER TOOL** – Used to move test samples from the bending clamp into the specimen holders.

To conduct a test, the material is cut to the required shape and shape and marked with a concentrated notch. An induced stress is made on the sample and it is place into the test sample rack and moved to the test tube which is filled with the test reagent and placed into the test bath at a constant temperature. At timed intervals, the samples are checked for any cracks developing at 90° to the notch. The number of failures is recorded.

## ENVIRONMENTAL STRESS CRACKING APPARATUS (RR/ESC)


### TECHNICAL SPECIFICATION

- Simple determination of Stress Cracking
- PID electronic temperature control
- Resolution 0.1°C
- PT100 PRT sensor accurate to 0.1°C
- Integrated stirrer motor
- Electronic timer HR:MIN:SEC
- Stainless steel liquid bath
- Sample cutter
- Sample nicking jig
- Sample bending clamp
- Sample specimen holders (Stainless Steel)
- 48 glass test tubes
- Rubber stoppers
- Product user manual
- CE declaration certificate
- 1 year return to base warranty
- Conforms to ASTM D1693

### OPTIONAL ANCILLARIES

- Sample specimen holders (Brass)

| WEIGHTS & DIMENSIONS: RR/ESC |    |
|------------------------------|----|
| Net Weight (kg)              | 25 |
| Width (cm)                   | 50 |
| Depth (cm)                   | 40 |
| Height (cm)                  | 50 |



## SAMPLE & COMPONENT TESTING

# UNRESTRAINED LINEAR THERMAL FILM SHRINKAGE LIQUID IMMERSION METHOD

Designed and manufactured by Ray-Ran, the **Unrestrained Linear Thermal Film Shrinkage** apparatus or more commonly referred to as the liquid immersion method, is used to determine the thermal shrinkage of plastic film & sheeting.

As a result of manufacturing processes, internal stresses may be locked into a film or sheet which can be released by heating which causes shrinkage of the material. The amount of shrinkage which takes place is dependent on the temperature during the test. The results are plotted on a graph showing percentage shrinkage against temperature.

Manufactured to a high standard, the apparatus is supplied with electronic digital temperature control, temperature bath and all tooling accessories. The integrated temperature bath fitted with digital temperature controller, heater and a PT100 platinum resistance thermometer accurately maintains the test temperature to within  $\pm 0.5$  °C. To ensure temperature



stability within the bath, a stirrer motor system is fitted, and the built in electronic timer monitors the duration of the test.

The apparatus comes as standard with stainless-steel specimen holders which hold the film samples of 100mm x 100mm square. A sample cutting template of the same size is also supplied to ensure that the test samples are accurately cut to the required test dimensions.

To conduct a simple test, the apparatus is heated to test temperature and the test samples immersed in the temperature bath medium. After 10 seconds they are removed and allowed to cool back to room temperature. A measurement of the sample is taken longitudinal and transverse directions and the percentage shrinkage calculated.

## THERMAL FILM SHRINKAGE LIQUID IMMERSION METHOD (RR/FSL)

### TECHNICAL SPECIFICATION

- Simple determination of Unrestrained Film Shrinkage
- Liquid immersion method
- Simple to operate
- PID electronic temperature control
- Resolution 0.1°C
- PT100 PRT sensor accurate to 0.1°C
- Integrated stirrer motor
- Electronic timer HR:MIN:SEC
- Stainless Steel liquid bath
- Sample cutting template 100mm x 100mm
- 6 off Sample specimen holders (Stainless Steel)
- Product user manual
- CE declaration certificate
- 1 year return to base warranty
- Conforms to ASTM D2732

### OPTIONAL ANCILLARIES

- Sample specimen holders (Stainless Steel)

| WEIGHTS & DIMENSIONS: RR/FSL |    |
|------------------------------|----|
| Net Weight (kg)              | 25 |
| Width (cm)                   | 50 |
| Depth (cm)                   | 40 |
| Height (cm)                  | 50 |



## SAMPLE & COMPONENT TESTING

# UNRESTRAINED LINEAR THERMAL FILM SHRINKAGE HOT PLATE METHOD

Designed and manufactured by Ray-Ran, the **Hot Plate Thermal Film Shrinkage** Apparatus is used to determine the thermal shrinkage of plastic film & sheeting and is ideal for most quality control procedures where precise results are essential.

The principal of the hot plate method to determine film shrinkage is similar to the liquid immersion method, notably to determine the amount of shrinkage of a film when heated to release the internal stresses produced due to the manufacturing methods.

The integrated hot plate with digital temperature controller, heater and a PT100 platinum resistance thermometer accurately maintains the test temperature to within  $\pm 0.1^\circ\text{C}$  and the built-in electronic timer monitors the duration of the test.

The hot plate test procedure is very simple and quick to conduct. A film sample of round or rectangular shape is placed on a hot plate which is maintained at a constant temperature. The surface of the hot plate is lightly oiled with silicone oil to ensure good heat transfer from the hot plate to the film and a light cover plate which keeps the film flat is placed on the top of the film sample. After a set time, the film is removed



from the hot plate and cooled to room temperature. It is then re-measured and its percentage shrinkage determined.

Although the Hot Plate Film Shrinkage Apparatus is not controlled by any international test standards, it's economical features and high temperature accuracy make it a valuable piece of testing equipment within the polymer industry.

## THERMAL FILM SHRINKAGE HOT PLATE METHOD (RR/FSHP)

### TECHNICAL SPECIFICATION

- Simple determination of Unrestrained Film Shrinkage
- Hot plate method
- Simple to operate
- PID electronic temperature control
- Resolution  $0.1^\circ\text{C}$
- PT100 PRT sensor accurate to  $0.1^\circ\text{C}$
- Electronic timer HR:MIN:SEC
- Sample cover plate
- Product user manual
- CE declaration certificate
- 1 year return to base warranty

### OPTIONAL ANCILLARIES

- Test sample cutter  $\varnothing 50\text{mm}$
- Test sample cutter  $50\text{mm} \times 50\text{mm}$
- Silicone oil (100ml bottle)

| WEIGHTS & DIMENSIONS: RR/FSHP |    |
|-------------------------------|----|
| Net Weight (kg)               | 6  |
| Width (cm)                    | 20 |
| Depth (cm)                    | 30 |
| Height (cm)                   | 25 |



# UNIVERSAL TESTING MACHINES X250

The **X250 Universal Testing Machine** is the smallest of the models offered by Ray-Ran within the product range. The small footprint single column machine has a range of 5N to 3KN and is used in all industries worldwide for materials testing, product testing, research and development, and quality control applications. Fully automated by advanced computer technology the testing system with high precision control ensures simplicity in operation and highly accurate test results.

The Universal Testing Machine uses feature-rich **WinTest™ Analysis** software for control and testing, with comprehensive onboard test methods that are compliant with all relevant industry standards including advanced, standard and complex multi-stage test procedures which are fully configurable and controlled using a standard PC interface. An optional IPC3 industrial-grade Panel PC with touchscreen control can be supplied making the machine a robust standalone system without the need for an external PC or laptop, or customers can utilize their own desktop computer or laptop to operate the machine.




84

## THE STANDARD FEATURES OF THE MACHINE INCLUDE:

- A fully digital testing system with high precision control and accuracy includes automated computer control of test methods giving simplicity of operation.
- High-resolution auto ranging load cells with accuracies better than +/-0.5% down to 1/1000th of the load cell capacity.
- Automatic recognition and calibration of load cells and extensometers, with instant calibration check facility.
- 800% overload capability of load cells without damage.
- Small footprint design, giving the economy of a bench and floor space.
- High efficiency pre-loaded self-cleaning ball screws for fast, quiet testing. Fitted with sealed for life lubricated end bearings.
- Crosshead guidance system providing precise alignment and smooth running.
- Precision crosshead control via digital AC servo drive and brushless servo motor giving maintenance free operation and 4,000,000 steps per revolution positional control.



- High-speed data collection systems for up to 4 synchronous channels.
- 6 I/O channels for additional devices such as extensometers, micrometers, calipers, balances, etc.
- High stiffness loading frames with solid specialized steel crossheads and rigid extruded support columns with T-slots for accessory mounting.
- Overload, over travel and impact protection.
- Telescopic covers giving additional protection for ball screws against dust and testing debris.
- Extensive range of grips and fixtures for tension, compression, flexural, shear, peel and product testing, etc.
- A wide range of contacting and non-contacting extensometers is available including laser and video models. Force Measurement Universally Calibrated, better than Grade 0.5 EN 7500-1, DIN 51221 ASTM E-4. AFNOR A03-501. Range 0.4% to 100% minimum.
- Automatic identification of load cell. Resolution 1 part in 500000 with auto-ranging. Electronic load cell protection.
- Extension Measurement. Full frame length to a maximum resolution 0.000001 (selectable). Accuracy ±0.01mm. Absolute, relative and auxiliary modes mm, in and percent. Programmable extension limits.



85

- Speed Control. Drive system temperature and current protection. Positional jog speed 0.00001mm/min.
- Load Frame. Rigid frame, using dual slide crosshead guidance system and rigid extruded support column. Frame stiffness 8kN/mm plus K factor facility built-in.
- Re-circulating ball screw with bellows.
- Electronic limit trips, total travel trips, and customer programmable safety stops.
- Rubber mat front protection.
- Remote connection self-diagnostics software and registration. Giving direct connection of your machine for remote service, on-line training, software upgrades, new test methods uploaded, etc


## TECHNICAL SPECIFICATION

- Single column machine bench mounted (optional base cabinet available)
- A/C servo drive
- Machine Capacity 5N – 3kN
- Auto Load cell identification
- Load cell protection
- High-resolution load cell
- User-selectable load and travel limits
- Speed Range 0.00001 to 2500mm/min in steps of 0.00001mm/min
- Crosshead Travel (excluding grips) 630-990mm
- Telescopic lead screw covers
- Throat 200mm
- Vertical space 800-1160mm
- Direct connection to extensometers
- **WinTest™ Analysis** operating software
- Ethernet (or USB via adaptor)
- Supply voltage 230V, 1ph 50/60Hz (115V option available)
- Power 0.3kW
- Product user manual
- CE declaration certificate
- 1-year return to base warranty

## OPTIONAL ANCILLARIES

- Range of manual grips
- Range of pneumatic grips with a footswitch
- Compression platens, round or square
- 3 & 4 point bending fixtures various spans available
- Up to 4 Load cells can be fitted to a maximum of 3KN
- Various Extensometers for different materials available
- Wide frame extension
- Height extension
- Remote connection software

| WEIGHTS & DIMENSIONS: RR/UTM X250 |     |
|-----------------------------------|-----|
| Net Weight (kg)                   | 65  |
| Width (cm)                        | 69  |
| Depth (cm)                        | 53  |
| Height (cm)                       | 149 |



# UNIVERSAL TESTING MACHINES X350

The **X350 Universal Testing Machine** is the middle range of the models offered by Ray-Ran within the product range. The small footprint double column machine has a range of 5N to 20KN and is used in all industries worldwide for materials testing, product testing, research and development, and quality control applications. Fully automated by advanced computer technology the testing system with high precision control ensures simplicity in operation and highly accurate test results.


The Universal Testing Machine uses feature-rich **WinTest™ Analysis** software for control and testing, with comprehensive onboard test methods that are compliant with all relevant industry standards including advanced, standard and complex multi-stage test procedures which are fully configurable and controlled using a standard PC interface. An optional IPC3 industrial-grade Panel PC with touchscreen control can be supplied making the machine a robust standalone system without the need for an external PC or laptop, or customers can utilize their own desktop computer or laptop to operate the machine.

## THE STANDARD FEATURES OF THE MACHINE INCLUDE:

- A fully digital testing system with high precision control and accuracy includes automated computer control of test methods giving simplicity of operation.
- High-resolution auto ranging load cells with accuracies better than +/-0.5% down to 1/1000th of the load cell capacity.
- Automatic recognition and calibration of load cells and extensometers, with instant calibration check facility.
- 800% overload capability of load cells without damage.
- Small footprint design, giving the economy of a bench and floor space.
- High efficiency pre-loaded self-cleaning ball screws for fast, quiet testing. Fitted with sealed for life lubricated end bearings.
- Crosshead guidance system providing precise alignment and smooth running.
- Precision crosshead control via digital AC servo drive and brushless servo motor giving maintenance free operation and 4,000,000 steps per revolution positional control.



- High-speed data collection systems for up to 4 synchronous channels.
- 6 I/O channels for additional devices such as extensometers, micrometers, calipers, balances, etc.
- High stiffness loading frames with solid specialized steel crossheads and rigid extruded support columns with T-slots for accessory mounting.
- Overload, over travel and impact protection.
- Telescopic covers giving additional protection for ball screws against dust and testing debris.
- Extensive range of grips and fixtures for tension, compression, flexural, shear, peel and product testing, etc.
- A wide range of contacting and on-contacting extensometers is available including laser and video models. Force Measurement Universally Calibrated, better than Grade 0.5 EN 7500-1, DIN 51221 ASTM E-4. AFNOR A03-501. Range 0.4% to 100% minimum.
- Automatic identification of load cell. Resolution 1 part in 500000 with auto-ranging. Electronic load cell protection.
- Extension Measurement. Full frame length to a maximum resolution 0.000001 (selectable). Accuracy ±0.01mm. Absolute, relative and auxiliary modes mm, in and percent. Programmable extension limits.
- Speed Control. Drive system temperature and current protection. Positional jog speed 0.00001mm/min.



- Load Frame. Rigid frame, using dual slide crosshead guidance system and rigid extruded support column. Frame stiffness 50kN/mm plus K factor facility built-in.
- Re-circulating ball screw with bellows.
- Electronic limit trips, total travel trips, and customer programmable safety stops.
- Rubber mat front protection.
- Remote connection self-diagnostics software and registration. Giving direct connection of your machine for remote service, on-line training, software upgrades, new test methods uploaded, etc

## TECHNICAL SPECIFICATION

- Dual column machine bench mounted (optional base cabinet available)
- A/C servo drive
- Machine Capacity 5N - 20kN
- Auto Load cell identification
- Load cell protection
- High-resolution load cell
- User-selectable load and travel limits
- Speed Range 0.00001 to 2000mm/min in steps of 0.00001mm/min
- Crosshead Travel (excluding grips) 1100mm
- Telescopic lead screw covers
- Throat 200mm
- Vertical space 1275mm
- Direct connection to extensometers
- **WinTest™ Analysis** operating software
- Ethernet (or USB via adaptor)
- Supply voltage 230V, 1ph 50/60Hz (115V option available)
- Power 0.45kW
- Product user manual
- CE declaration certificate
- 1-year return to base warranty

## OPTIONAL ANCILLARIES

- Range of manual grips
- Range of pneumatic grips with a footswitch
- Compression platens, round or square
- 3 & 4 point bending fixtures various spans available
- Up to 4 Load cells can be fitted to a maximum of 20KN
- Various Extensometers for different materials available
- Wide frame extension
- Height extension
- Remote connection software

| WEIGHTS & DIMENSIONS: RR/UTM X350 |     |
|-----------------------------------|-----|
| Net Weight (kg)                   | 120 |
| Width (cm)                        | 69  |
| Depth (cm)                        | 53  |
| Height (cm)                       | 158 |



# UNIVERSAL TESTING MACHINES X500



The **X500 Universal Testing Machine** is the largest of the models offered by RayRan within the product range. The small footprint double column machine has a range of 5N to 100kN and is used in all industries worldwide for materials testing, product testing, research and development, and quality control applications. Fully automated by advanced computer technology the testing system with high precision control ensures simplicity in operation and highly accurate test results.

The Universal Testing Machine uses feature-rich **WinTest™ Analysis** software for control and testing, with comprehensive onboard test methods that are compliant with all relevant industry standards including advanced, standard and complex multi-stage test procedures which are fully configurable and controlled using a standard PC interface. An optional IPC3 industrial-grade Panel PC with touchscreen control can be supplied making the machine a robust standalone system without the need for an external PC or laptop, or customers can utilize their own desktop computer or laptop to operate the machine.

## THE STANDARD FEATURES OF THE MACHINE INCLUDE:

- A fully digital testing system with high precision control and accuracy includes automated computer control of test methods giving simplicity of operation.
- High-resolution auto ranging load cells with accuracies better than +/-0.5% down to 1/1000th of the load cell capacity.
- Automatic recognition and calibration of load cells and extensometers, with instant calibration check facility.
- 800% overload capability of load cells without damage.
- Small footprint design, giving the economy of a bench and floor space.
- High efficiency pre-loaded self-cleaning ball screws for fast, quiet testing. Fitted with sealed for life lubricated end bearings.
- Crosshead guidance system providing precise alignment and smooth running.
- Precision crosshead control via digital AC servo drive and brushless servo motor giving maintenance free operation and 4,000,000 steps per revolution positional control.
- High-speed data collection systems for up to 4 synchronous channels.



- 6 I/O channels for additional devices such as extensometers, micrometers, calipers, balances, etc.
- High stiffness loading frames with solid specialized steel crossheads and rigid extruded support columns with T-slots for accessory mounting.
- Overload, over travel and impact protection.
- Telescopic covers giving additional protection for ball screws against dust and testing debris.
- Extensive range of grips and fixtures for tension, compression, flexural, shear, peel and product testing, etc.
- A wide range of contacting and non-contacting extensometers is available including laser and video models. Force Measurement Universally Calibrated, better than Grade 0.5 EN 7500-1, DIN 51221 ASTM E-4. AFNOR A03-501. Range 0.4% to 100% minimum.
- Automatic identification of load cell. Resolution 1 part in 500000 with auto-ranging. Electronic load cell protection.
- Extension Measurement. Full frame length to a maximum resolution 0.000001 (selectable). Accuracy ±0.01mm. Absolute, relative and auxiliary modes mm, in and percent. Programmable extension limits.
- Speed Control. Drive system temperature and current protection. Positional jog speed 0.00001mm/min.
- Load Frame. Rigid frame, using dual slide crosshead guidance system and rigid extruded support column. Frame stiffness up to 400kN/mm plus K factor facility built-in.
- Re-circulating ball screw with bellows.
- Electronic limit trips, total travel trips, and customer programmable safety stops.
- Rubber mat front protection.
- Remote connection self-diagnostics software and registration. Giving direct connection of your machine for remote service, on-line training, software upgrades, new test methods uploaded, etc

## TECHNICAL SPECIFICATION

- Dual column machine bench mounted (optional base cabinet available)
- A/C servo drive
- Machine Capacity 5N - 100kN
- Auto Load cell identification
- Load cell protection
- High-resolution load cell
- User-selectable load and travel limits
- Speed Range 0.00001 to 2000mm/min in steps of 0.00001mm/min
- Crosshead Travel (excluding grips) 1100mm
- Telescopic lead screw covers
- Throat 200mm
- Vertical space 1300mm
- Direct connection to extensometers
- **WinTest™ Analysis** operating software
- Ethernet (or USB via adaptor)
- Supply voltage 230V, 1ph 50/60Hz (115V option available)
- Power 1.2kW
- Product user manual
- CE declaration certificate
- 1-year return to base warranty

## OPTIONAL ANCILLARIES

- Range of manual grips
- Range of pneumatic grips with a footswitch
- Compression platens, round or square
- 3 & 4 point bending fixtures various spans available
- Up to 4 Load cells can be fitted to maximum of 100kN
- Various Extensometers for different materials available
- Wide frame extension
- Height extension
- Remote connection software

| WEIGHTS & DIMENSIONS: RR/UTM X500 |     |
|-----------------------------------|-----|
| Net Weight (kg)                   | 235 |
| Width (cm)                        | 78  |
| Depth (cm)                        | 53  |
| Height (cm)                       | 177 |

# WINTEST™ ANALYSIS

## THE ULTIMATE SOFTWARE PACKAGE IN FORCE MEASUREMENT

**WinTest™ Analysis** is the advanced software suite supplied with the Ray-Ran range of **Universal Testing Machines**. It is a multi-functional and fully customisable software package that supports all industry standards including ISO, ASTM and BS EN specifications. The software covers tensile, compression, peel, shear, tear, cyclic, creep and multistage testing. It includes a wide range of industry standard test methods and the facility to create and store an unlimited number of bespoke test methods. Additional flexibility is provided by user-defined multistage step testing for highly specialised testing requirements and includes an extensive range of calculations applicable to many industries, including all variations of force, elongation, stress and strain values and many others.

The virtual control panel allows the operator full control of all tester functions and the ability to conduct simple tests manually as well as providing easy access to stored test methods, system configuration and diagnostics. Test results are displayed in real time on the auto scaling graph and in tabular format with each test in the batch highlighted in a different colour. The software also offers automated storage of all test data and test results.



Test reports are easily generated, and test data can be transferred directly to other software applications for enhanced report generation, SPC trend analysis or laboratory manager reports for example. They can be exported to Microsoft Word and/or Excel to provide you with full editing features and copy and paste capability to produce presentation-quality test reports, charts or test data in spreadsheet format. Including company details and logos to complete that professional look.

Test reports are easily converted into Adobe PDF files so you can simply email them as a PDF attachment or

send them as a Microsoft Word document or an Excel file. Test data can also be exported in ASCII delimited format to Microsoft Excel or other spreadsheet software using search criteria such as References, Start Date and End Date. It is also possible to export Minimum, Mean and Maximum results data based on a periodic selection (daily, weekly, monthly etc.) to highlight trends over a specified interval.

The integrated HTML help file with added search function includes simple explanations of machine operation, test result descriptions and graphical Flash

representation of tests and test calculations. This feature also enables you to view graphically how specific test results are calculated to help you verify the correct selection of test calculations.

### RESULTS FEATURES INCLUDE:

- Comprehensive library of industry standard calculations.
- Configurable statistics summary for each test report.
- Customised test calculations.
- Visual Pass/Fail.
- Sequential calculations to take measurement at set intervals for long term tests etc.
- Custom statistics can be generated for selected calculation.
- Retrospective analysis of all test calculations.
- Calculate results on pre-defined test regions.
- Comments field and custom columns available for each test series and for individual tests.
- User-friendly test data backup can be configured for periodic reminders.
- Free backup / archive viewer program available

### REPORT FEATURES INCLUDE:

- Fully configurable test reports
- User-defined header and footer on test reports.

- Generate test reports in PDF format for email etc.
- Auto-print and preview option.
- Report transfer into Microsoft Excel™ and Microsoft Word™.
- Import and export of test definitions in XML format
- Export test results in ASCII format
- Export raw curve data

### GRAPH FEATURES INCLUDE:

- Pass/Fail tolerance bands
- Display of best fit straight line in the elastic region, for calculation of E modulus, proof stress etc.
- Golden sample, a test curve can be selected as a reference and tolerance bands can be set to provide an instant visual check that all subsequent tests are within tolerance.
- Event marking during real time plotting of test curve.
- Visual display of calculations.
- User defined annotations.
- Built-in drawing function.
- Visual event markers.
- User defined peaks and troughs.
- Show the current graph co-ordinates for the current mouse position

## WINTEST™ REPORTS

**WinTest™ Reports** is an enhancement to WinTest Analysis to add flexibility to data analysis and statistical reporting. The package provides a report generation capability that can include long-term statistics and control charts for all specified calculations. **WinTest™ Reports** can also be configured to display headings, titles, company logos, graphs, charts, pop-up menus and specific technical information.



SAMPLE & COMPONENT TESTING



DM-T DIGITAL MICROMETER

The **DM-T Digital Micrometer** available from Ray-Ran has been specifically designed to quickly measure with precision and accuracy, the thickness of a variety of materials including foils, paper, natural or synthetic fabrics, leathers, metals, plastics and rubber. Operation is simple and the large LCD display shows the thickness of the materials under test to a resolution of 0.5 microns and an accuracy of 1.0 micron.



The motor-driven instrument utilizes the dead weight micrometer principle for high accuracy and repeatability and has a lower anvil and movable pressure foot which is made from lapped, stainless steel. Special anvil diameters and dead weight loads are also available.

The Apparatus has user configurable options for altering speed up to 2.5mm/sec, dwell time 100 msec – 6 sec, time of day, zero threshold, live readings, starting units, LED pulse, passwords and calibration menus. The apparatus is fitted with an RS232 serial connector so it can be used with the optional GraphMaster™ software for results presentation and storage.

**The instrument can be configured to meet any of the standards listed below:**  
TAPPIT 411, ISO 534, DIN 53105, ISO 3034, APPITA 1301.426, PAPTAC D.4, ISO 037, ISO 4593(\*), ASTM D6988(\*), ASTM F2251, ISO 5084

DIGITAL MICROMETER (RR/DM-T)

TECHNICAL SPECIFICATION

- Resolution: 0.5 microns
- Accuracy: 1.0 microns
- Measurement Range: 0 – 0.05” (0 – 1.27mm)
- 0 – 0.50” (12.7mm) available for thicker samples.
- Output: RS232
- Power: 110V/220V 50Hz/60Hz
- User configurable for units, dwell time and down speed and more
- Measurement units in microns, mils, mm, inches
- Easy to use touch button
- Lapped stainless steel pressure foot and anvil
- Product user manual
- CE declaration certificate
- 1 year return to base warranty

OPTIONAL ANCILLARIES

- GraphMaster™ Software
- Special anvil diameters and dead weight loads
- ¼” dia anvil (minimum pressure 8 PSI)
- ½” dia anvil (minimum pressure 4 PSI)

| WEIGHTS & DIMENSIONS: MODEL RR/DM-T |      |
|-------------------------------------|------|
| Net Weight (kg)                     | 10   |
| Width (cm)                          | 15   |
| Depth (cm)                          | 24   |
| Height (cm)                         | 30.5 |

SAMPLE & COMPONENT TESTING



DM-B DIGITAL MICROMETER

The **DM-B Digital Micrometer** available from Ray-Ran combines unmatched accuracy and resolution with a modern contemporary look. The Micrometer can be configured to meet ISO, ASTM, TAPPI, EDANA or other international specifications. The Micrometer offers a cantilever balance system to allow low pressure measurements. This feature also allows adding or removing additional weights for multiple pressure applications. The instrument is supplied with a sensing eye next to the anvil so that when a sample is detected the test cycle automatically starts.

The rigid frame and linear encoder enable repeatable and accurate results. A large touchscreen displays test data to 0.1µm resolution. Individual and rolling average results are displayed after each test. A statistic mode provides further analysis. The software with intuitive menu allows easy operation of the apparatus.

Testing is quick and easy. Parameters are software controlled and can be adapted, such as gap height, auto test mode, auto-zero and dwell time. A cantilever measuring mechanism allows



for very light pressures so compressible materials such as thin films and paper tissue can be measured accurately without damaging the material under test. A variety of pressure foot weights are available to increase the pressure on the upper anvil and an optional combination anvil allows the operator to change the diameter of the pressure foot so different standards can be tested on one unit.

For ease of operation the apparatus can be used with an optional strip feeder and support bench for continuous testing and an optional foot switch.

The apparatus is fitted with an RS232 serial connector so it can be used with the optional GraphMaster™ software for results presentation, storage or file transfer directly to Excel®, Access®, or clipboard where test reports are automatically generated.

Applications include Paper, corrugated, cloth, plastic, plastic film, textile fabrics, nonwovens, battery separators, felts, leather, tissue paper and others.

**The apparatus conforms to**  
ASTM D374, D1777, D5729, D6988, ISO 534, 3034, 4593:1993, 5084, 9073-2, 12625-3, APPITA 1301.426, TAPPI T-411, EDANA 30.4, PAPTAC D.4, DIN 53370, BS 2782-6, WSP 120.1, WSP 12

DIGITAL MICROMETER (RR/DM-B)

TECHNICAL SPECIFICATION

- Resolution: 1.0/0.1 µm (0.01 mil/0.001 mil) on request
- Accuracy: within 0.001 mm (0.00004 in.) or 1% of paper thickness, Which ever is better
- Measurement Range: 0 - 10 mm; 0 - 10,000 µm
- Lowering Speed: 0.8 -5.9 mm/sec
- Output: RS232
- Power: 90V/230V 50Hz/60Hz
- User configurable for units, dwell time and down speed and more
- Measurement units in microns, mils, mm, inches
- Easy to use touch screen
- Selectable pressure foot and anvils for different pressures
- Product user manual
- CE declaration certificate
- 1 year return to base warranty

OPTIONAL ANCILLARIES

- GraphMaster™ Software
- Special anvil diameters and pressures
- Strip feeder
- Support plate
- Foot switch

| WEIGHTS & DIMENSIONS: MODEL RR/DM-B |      |
|-------------------------------------|------|
| Net Weight (kg)                     | 10   |
| Width (cm)                          | 15   |
| Depth (cm)                          | 24   |
| Height (cm)                         | 30.5 |

## SAMPLE & COMPONENT TESTING

# INDENTATION HARDNESS APPARATUS

For accurately checking the shore hardness of a wide range of elastomers and plastic materials, the Ray-Ran range of **Durometers** are simply hard to beat. Since material hardness is a key characteristic of an elastomer and one which can be related to important material properties, it is beneficial to be able to accurately measure the material hardness of a supplied product making them ideal for product design and quality control inspection.

Offered in analogue or digital scales, each durometer can be used handheld or fitted to the Ray-Ran bench mounted Lever Operated Durometer Stand. Each durometer available utilizes the Shore method of hardness measurement and features a spring-loaded indenter module which



is pressed down squarely on the surface of the material to be measured resulting in its shore hardness.

The unique interchangeable modules allow for a wide and useful range of the Shore harnesses scales to be measured on one instrument. Depending on the material being tested, hardness scales are available to Shore A, B, C, D, DO and OO with each interchangeable module supplied calibrated to the durometer head. The only exception to this interchangeability is the softest Shore OO range due to its very low indenter loads.

The models offered within the Ray-Ran range are the STD226 Digital Durometer and the STD227 Analogue Durometer. Each durometer comes with its own calibration certificate tested to the relevant module, protective carry case, product user manual and 12 months return to base warranty.

## INDENTATION HARDNESS (RR/STD226, RR/STD227)

### TECHNICAL SPECIFICATION

#### MODEL STD 226 DIGITAL DUROMETER

- Resolution - half of one degree Shore hardness scale
- Easy to use - hand held
- Sturdy, accurate and reliable
- BCD output facility
- Uses 1 x SR44 flat battery
- Multiple Modules with the same instrument head
- Supplied in convenient storage/protective case

#### MODEL STD 227 ANALOGUE DUROMETER

- Resolution - one degree Shore hardness scale
- Supplied complete with tolerance settings
- Optional maximum hand
- Multiple Modules with the same instrument head
- Supplied in convenient storage/protective case

### OPTIONAL ANCILLARIES

- Lever operated bench stand
- Set of 6 coloured rubber reference blocks (shore A only)
- Shore A Scale ASTM D 2240, ISO 7619 & 868, DIN 53505
- Shore B Scale ASTM D 2240
- Shore C Scale ASTM D 2240
- Shore D Scale ASTM D 2240, ISO 7619 & 868, DIN 53505
- Shore DO Scale ASTM D 2240
- Shore O Scale ASTM D 2240

| WEIGHTS & DIMENSIONS:<br>INDENTATION HARDNESS APPARATUS |        |        |
|---|--------|--------|
|   | STD226 | STD227 |
| Net Weight (kg)   | 0.3    | 0.22   |
| Width (cm)  | 6.5    | 3.8    |
| Depth (cm)  | 2.5    | 2.3    |
| Height (cm)   | 15     | 13.5   |

## LEVER OPERATED DUROMETER BENCH STAND (RR/DBS)

### TECHNICAL SPECIFICATION

- Can be fitted with STD227 and STD226 durometers
- Easy to use lever operated
- Large sample lever platform
- Adjustable durometer height
- Set of test weights supplied
- "V" block for cylindrical samples
- Product user manual
- 12 month return to base warranty

### OPTIONAL ANCILLARIES

- STD226 Durometer
- STD227 Durometer

| WEIGHTS & DIMENSIONS: RR/DBS |    |
|------------------------------|----|
| Net Weight (kg)              | 5  |
| Width (cm)                   | 35 |
| Depth (cm)                   | 25 |
| Height (cm)                  | 45 |

## SAMPLE & COMPONENT TESTING

# LEVER OPERATED DUROMETER BENCH STAND

The **Lever Operated Durometer Bench Stand** has been designed by Ray-Ran to enable the simple mounting of a Durometer to be brought into contact with a test sample in a vertical direction with a known constant loading force being applied. This almost eliminates any operator influence on the results, thus ensuring greater accuracy and repeatability.

Simple to operate, the durometer is loaded with a specified static load and the test sample is placed on the lever operated surface plate. The test sample is then pressed against the durometer and the constant force of the static load applied ensuring accurate repeatable test results.

Manufactured to comply with international testing methods ASTM D2240 and DIN 53505, the apparatus can accommodate both the STD226 and STD227 durometers. The bench stand is supplied complete with all test weights to the relevant test standards and "V" block for testing polymer tube samples.



## DEMONSTRATION FACILITY

Ray-Ran Test Equipment is very pleased to welcome customers to its purpose-built onsite demonstration facility based in Nuneaton. The facility is ideal for product verification and validation purposes which lets customers get a “hands on” experience of the equipment of interest to ensure their future investment meets their requirements.

For customers who have already purchased equipment from Ray-Ran the onsite demonstration facility also provides an excellent place to train your personnel. Our fully trained technicians ensure your operators are fully conversant with all the technical aspects of the equipment and understand the operation procedures ensuring accurate test results maximizing the return on your investment.

Based within a short distance of Birmingham Airport, Train and Motorway routes, Ray-Ran Test Equipment is very accessible to many UK and European customers. If your journey or visit requires an overnight stay, Ray-Ran can organise local hotel accommodation if required.

If you feel this facility would be beneficial to you, please contact Ray-Ran to arrange your visit.

WE LOOK FORWARD TO WELCOMING YOU TO RAY-RAN TEST EQUIPMENT LTD.

If an onsite demonstration is not possible simply send us your test samples and we will conduct a small batch test on the equipment of interest to verify the apparatus meets your testing requirements.

The following products are always available for demonstration. For other products, please check in advance as they may be subject to a lead time.

- Pendulum Impact and Puncture Tester
- Melt Flow Indexer
- Pneumatic and Hand Cutting Presses
- CNC1 Milling Machines
- HDT/Vicat
- Test Sample Moulding Apparatus
- Density Balance
- Density Gradient Colum
- Apparent Bulk Density
- Hardness Testers
- Coefficient of Friction Tester
- Sample Notching Cutter
- Falling Dart

## SERVICE & CALIBRATION

Ray-Ran recommend a regular service and calibration programme for your testing equipment after purchase to maintain their accuracy and reliability. All service and calibration procedures are conducted by skilled technicians with the most up to date equipment available ensuring your apparatus fully complies with International Testing Standards. Calibrations are usually conducted on a 24hr turnaround basis ensuring limited downtime and expense to the customer. Each machine is supplied with a Certificate of Calibration which is your guarantee of our competence and of the continuing accuracy of your instrumentation.

## REPAIRS

In the event that your equipment breaks down or develops a fault Ray-Ran offer a complete repair program. Conducted by fully skilled technicians the work will include fault diagnosis followed by a plan of action to make the equipment fully functional. Ray-Ran will notify you of all the costs involved prior to work being carried out and all repairs are covered with a 12-month warranty.

## COMMITMENT TO QUALITY

Ray-Ran Test Equipment Ltd is an ISO 9001:2015 UKAS accredited company and is proud of its record in supplying high quality products to its customers. The Company’s policy of commitment and dedication to support its customers ensures Ray-Ran remains the world’s leader in supplying polymer and materials testing equipment.

# TESTING INSTRUMENTS

HIGH PRECISION TESTING INSTRUMENTATION

# SERVICE AND CALIBRATION

WORLDWIDE INSTALLATION, MAINTENANCE  
AND CALIBRATION SERVICE

# TECHNICAL ASSISTANCE

EXPERT ENGINEERING AND APPLICATIONS SUPPORT

# PRODUCT DESIGN & DEVELOPMENT

EXPERT ON SITE R & D, BESPOKE AND 1 OFF CAPABILITY.  
CONTACT US WITH YOUR REQUIREMENTS.

## LAB SOLUTIONS

9 Cours Manuel de Falla  
Z.I. MOZART  
26000 Valence  
FRANCE  
+33 04 75 60 11 77  
[www.labsolutions.fr](http://www.labsolutions.fr)

AGENT



## RAY-RAN TEST EQUIPMENT LTD

KELSEY CLOSE,  
ATTLEBOROUGH FIELDS  
INDUSTRIAL ESTATE  
NUNEATON,  
WARWICKSHIRE  
CV11 6RS

TEL: +44 (0)24 7634 2002  
FAX: +44 (0)24 7664 1670  
EMAIL: [polytest@ray-ran.com](mailto:polytest@ray-ran.com)  
WEB: [www.ray-ran.com](http://www.ray-ran.com)

# Ray / Ran